



# AVENIO Oncology Analysis Software

Exceptional performance, minimal effort



# AVENIO Oncology Analysis Software


The **AVENIO Oncology Analysis Software** processes sequencing data from plasma and tissue samples run using the AVENIO ctDNA Analysis Kits V2 and AVENIO Tumor Tissue Analysis Kits V2.

The interface guides users confidently through the post-sequencing workflow. Intelligent bioinformatics with advanced proprietary algorithms and error suppression strategies deliver comprehensive results with proven accuracy and reproducibility.<sup>1,2</sup>


Powerful and intuitive, the AVENIO Oncology Analysis Software offers simplicity, efficiency, and accuracy.




Log-in




Select plasma or tissue




Run settings



Analysis



Results




Reporting

## Set up and run the analysis with ease



The AVENIO Oncology Analysis Software requires minimal expertise to set up and run an analysis. This allows laboratory and bioinformatics professionals to focus more of their time on the findings and insights, and less on managing information systems.

Intuitive user interface:

- Select plasma or tissue
- Offers two user types (Admin, Lab User)
- Provides real-time tracking of analysis status
- Enables filtering, custom annotation, and analytical concordance reporting between tissue and plasma
- Provides reports in PDF format, as well as variant results in VCF and BAM file formats.



With its ability to process approximately 800 million reads in less than a day, the AVENIO Oncology Analysis Software can be set up to run in the evening, delivering results ready for review in the morning.



## Detect all four mutation classes in a single DNA workflow

The software's bioinformatic algorithms and checkpoints have been optimized to enable accurate variant calls across all four mutation classes.<sup>1</sup>

## Detect low frequency alleles with exceptional sensitivity and specificity

The AVENIO Oncology Analysis Software employs integrated digital error suppression (iDES) strategies in the ctDNA workflow, combining molecular barcodes with in silico error suppression techniques. This proprietary method enables detection of low frequency alleles down to 0.1% with exceptional sensitivity and specificity.<sup>1,2</sup>

The tumor tissue workflow employs enzymatic error suppression and molecular barcodes for accurate detection of all four mutation classes from DNA only.

**AVENIO Oncology Analysis Software**  
Variant Report

**Sample Details**

SAMPLE ID	SNV1-flag2	ANALYSIS COMPLETION DATE	19 Apr 2018 03:22:00 PDT
SAMPLE TYPE	Plasma	PANEL	AVENIO ctDNA Targeted Panel
SAMPLE ADAPTER	SA 7	FILTER SET	Unfiltered
CUSTOM ANNOTATION	Plasma, ctDNA-avenio-special	INPUT DNA MASS	40.00 ng
PLASMA VOLUME	4.4 ml	ISOLATED DNA MASS	40.00 ng
SUBJECT ID	rowd.apc01at [REDACTED]	DATE OF SAMPLE COLLECTION	22 Jun 2017

**Results Summary**

Gene	Variant	Variant Description	Allele Fraction	No. of Mutant Molecules per mL
ALK	p.Tyr1827Cys	Missense variant	7.94%	262
ALK	p.Tyr287Cys	Missense variant	1.37%	45.2
ALK	c.3939-2A>C	Splice acceptor variant & Intron variant	1.37%	45.2
ALK	c.819-2A>C	Splice acceptor variant & Intron variant	1.37%	45.2
ALK	p.Gln1125Gly	Synonymous variant	8.84%	292
ALK	p.Gly85Gly	Synonymous variant	0.05%	1.60
ALK	p.Ser1081Cys	Missense variant	80.02%	2640
ALK	p.Ser110Cys	Missense variant	1.34%	44.2
ALK	p.Gly845Gly	Synonymous variant	0.21%	7.08
ALK	p.Asn767Ser	Missense variant	0.28%	9.14
ALK	p.Gly692Trp	Missense variant	86.54%	2860
ALK	p.Ser991Arg	Missense variant	28.95%	960
ALK	p.Gln500Gln	Synonymous variant	1.85%	60.9
ALK	p.Val476Ala	Missense variant	1.55%	50.5
APC	c.-2G>T	5 prime UTR premature start codon gain variant	0.98%	32.2
APC	c.-2G>T	5 prime UTR variant	24.11%	795
APC	p.Leu172Arg	Missense variant	31.96%	1060
APC	c.1912-5G>A	Splice region variant & Intron variant	2.30%	76.1
APC	p.Tyr486Tyr	Synonymous variant	0.10%	3.14
APC	p.Ala545Ala	Synonymous variant	0.14%	4.67
APC	p.Arg653Arg	Splice region variant & Synonymous variant	0.13%	4.26
APC	p.Ser837*	Stop gained	8.68%	287
APC	p.Tyr1376Asn	Missense variant	33.38%	1100
APC	p.Lys1454Glu	Missense variant		
APC	p.Thr1493Thr	Synonymous variant		

Sample ID: SNV1-flag2  
Report generation date: 19 Apr 2018 03:22:55 PDT  
FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

Page 1 of 161

# Performance metrics<sup>1</sup>

## AVENIO ctDNA Analysis Kits

Mutant Allele Frequency/ Copy Number	SNVs		Indels		Fusions		CNVs**
	0.5%*		1.0%*		1.0%		
AVENIO Tissue Kit	Sensitivity	PPV	Sensitivity	PPV	Sensitivity	PPV	Sample Dependent
Targeted	> 99%	> 99%	> 99%	> 99%	> 99%	> 99%	
Expanded	> 99%	> 98%	> 99%	> 99%	> 99%	> 99%	
Surveillance	> 99%	> 99%	> 99%	> 99%	> 99%	> 99%	

\*The reported sensitivity and PPV is for 50ng DNA input at the stated AF%. The kit has the ability to report variants down to 0.1% AF (with lower sensitivity), and can make libraries with as low as 10 ng.



\*\*The kit reports ERBB2, EGFR, and MET CNVs. Performance is dependent on various factors, including sample type, input DNA amount, percent tumor content, and the individual tumor's genomic copy number. Shown are examples of detected CNVs at ddPCR verified copy numbers, using a commercially available reference sample (50 ng). CNV specificity was high (PPV >99%). ERBB2 detected at 2.2 copies, MET detected at 2.6 copies.

Performance samples: Seraseq ctDNA Complete Mutation Mixes (SeraCare), healthy donor cfDNA. Sensitivity and Positive Predictive Value (PPV) metrics based on observed product performance. Sensitivity and PPV performance reported per variant. Sensitivity was determined using commercially available reference samples containing verified mutations at the stated allele frequencies. SNV performance data based on hotspot calls. The AVENIO ctDNA Analysis Kit V2 can achieve >99.99% per base specificity across each of the panels. Stated performance requires at least 40 million reads per sample for all AVENIO ctDNA Panels V2. Sequencing performed on an Illumina NextSeq 500/550/550Dx (Research Use Only mode) instrument

## AVENIO Tumor Tissue Analysis Kits

Mutant Allele Frequency/ Copy Number	SNVs		Indels		Fusions		CNVs
	5.0% *		5.0%		5.0%		
AVENIO ctDNA Kit	Sensitivity	PPV	Sensitivity	PPV	Sensitivity	PPV	Sample Dependent
Targeted	> 99%	> 99%	> 99%	> 99%	> 99%	> 99%	
Expanded	> 99%	> 98%	> 99%	> 99%	> 99%	> 99%	
Surveillance	> 99%	> 99%	> 99%	> 99%	> 99%	> 99%	

Samples: FFPE tissue curls/sections, DNA input: ≥20 ng of FFPE DNA, total DNA amount for each sample determined by input QC. Sensitivity and Positive Predictive Value (PPV) metrics based on observed product performance. Sensitivity and PPV performance reported per variant. Sensitivity was determined using commercially available reference samples containing verified mutations at the stated allele frequencies. SNV performance is panel wide. Indel, Fusion based on whitelist variants. AVENIO Tumor Tissue Analysis Kits V2 achieve >99.99% per base specificity across each of the panels. Stated performance requires at least 20 million reads per sample for Targeted, Expanded and Surveillance Kits V2. Sequencing performed on an Illumina NextSeq 500/550/550Dx (Research Use Only mode) instrument. ERBB2 detected at 4.5 copies, MET detected at 4.5 copies.



## Call mutations with confidence

Five leading oncology databases, a curated loci of interest list and a customizable annotation database are integrated into the AVENIO Oncology Analysis Software. This allows quick and easy access to trusted public resources for cross-referencing and verification of results:

- COSMIC
- TCGA
- ExAC
- dbSNP
- 1000 Genomes

**These annotation databases can help users call mutations with increased confidence.**

Published by:

**Roche Sequencing Solutions, Inc.**  
4300 Hacienda Drive  
Pleasanton, CA 94588

For Research Use Only. Not for use in diagnostic procedures.

1. Data on file with Roche.  
2. Newman AM, Lovejoy AF, Klass DM, et al. Integrated digital error suppression for improved detection of circulating tumor DNA. Nature Biotechnology. 2016;34(5):547-555. doi:10.1038/nbt.3520.

AVENIO is a trademark of Roche. NEXTSEQ is a trademark of Illumina. All other product names and trademarks are the property of their respective owners.

© 2025 Roche Sequencing Solutions, Inc. All rights reserved.

## Control quality across the workflow

The AVENIO Oncology Analysis Software provides a sequencing quality report with key metrics such as sequencing depth, number of reads, on-target rate, coverage uniformity, and error rate to confirm quality at different stages in the workflow.