

Roche



Roche Digital LightCycler® System

A technological guide to the powerful new addition to our PCR ecosystem

Digital 
LightCycler®

It's time for a leap forward in digital PCR technology.

Experience sensitivity, precision, flexibility, and integration in one powerful clinical research tool. The Digital LightCycler® System from Roche is a digital PCR system that can help laboratories push forward the boundaries of clinical research and has the potential to advance global medical knowledge.



Partitioning Engine

Height 9.8 in **Width** 9.8 in **Depth** 11.8 in

Touch-screen operation

Stand-alone to accommodate multi-room configuration

Partitioning Fluid

Inert non-volatile silicone fluid increases reliability and replicability and minimizes the chance of amplicon contamination.

3 unique nanowell plates

Standard SBS/96 MWP format

Height 128mm **Width** 85mm

8 sample lanes per plate



20,000 partitions High Sensitivity

166x87x160µm, ~45µL

Cell-free DNA

Oncology

Rare Mutation Detection



28,000 partitions Universal

121x62x128µm, ~30µL

Gene Expression

Transplant Rejection




100,000 partitions High Resolution

54x27x75µm, ~15µL

Copy Number Variation

NIPT

Human Genetic Disease



Analyzer

Height 35.4 in **Width** 39.4 in **Depth** 23.6 in

Fully integrated thermal cycling & partitioning imaging

Image capture within sealed nanowell plate

6 optical channel detection for multiplexing capabilities


Flexible batch size (increments of 8 up to 96 samples per run)

12-plate capacity

Choose the Digital LightCycler® System and unleash the true power of of digital PCR.

Run times*

Designed to use less overall time when there are fewer plates in the analyzer, the Digital LightCycler® System is also capable of processing large numbers of samples in an extremely fast average time per sample.



Volumes†

With average lost volumes of 10% on the 20,000 and 28,000 partition plates and just 5% on the 100,000 plate during research, the maximum waste volume is also extremely low.



Compatible optical dyes



Performance data

Quantification Precision* ≤5% for optimal sample input and ≤10% for low sample input

Quantification Accuracy +/-10% to the reference standard for optimal sample input and +/-20% for low sample input

Linearity and Dynamic Range At least 4-log of linear range with deviation from linear fit <0.2 on a log scale

CNV Assay Performance Discriminate 10% difference in CN on High Resolution Plate (100,000)

Rare Mutation Assay Performance LoD of 0.1% MAF on Universal Plate (28,000)

Indel Assay Performance LoD of 0.2% MAF on High Sensitivity Plate (20,000)

* Roche data on file: DH_02365.01_031B_Digital_LightCycler_Reagent_Feasibility_Report_v3, Document Number: 000000000001200000501942

Roche data on file: DH-02365.01-500E_Digital_LightCycler_System_Performance_report

Roche data on file: DHF_Digital_LightCycler_Reagent_Feasibility_Report

Roche data on file: DH-02365.01-008B_Digital_LightCycler_System_Feasibility_Report

† Assuming the Primer/Probe is at 5x (usually higher at 10x and 20x), master mix at 5x, restriction enzyme volume very small (neglected here)

* Quantified by the coefficients of variation (CV) of technical replicates



sequencing.roche.com

Published by:

Roche Sequencing and Life Science

915 Hague Road

Indianapolis, IN 46256 USA

<https://go.roche.com/dpcr>

Digital LightCycler® is a trademark of Roche.

© 2022 Roche Sequencing and Life Science, Inc. All rights reserved

MC-US-11707 BRO705001 A728 02/23