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1. Sano T, Oyama T, Kashiwabara K, et al. Expression status of p16 protein is associated with human papillomavirus oncogenic potential in cervical and genital lesions. *Am J Pathol.* 1998;153(6):1741-1748.
2. Negri G, Egarter-Vigl E, Kasal A, et al. p16INK4a is a useful marker for the diagnosis of adenocarcinoma of the cervix uteri and its precursors: a tissue microarray study with immunocytochemical correlations. *Am J Surg Pathol.* 2003;27(2):187-193.
3. von Knebel Doeberitz M, Vinokurova S. Host factors in HPV-related carcinogenesis: cellular mechanisms controlling HPV infections. *Arch Med Res.* 2009;40(6):435-442.
4. Vootoeve PM, Agarni R. The tumor-suppressive functions of the human INK4a locus. *Cancer Cell.* 2003;4:311-9.
5. Klaes R, Friedrich T, Spitkovsky D, et al. Overexpression of p16 (INK4a) as a specific marker for dysplasia and neoplastic epithelial cells of the cervix uteri. *Int J Cancer.* 2001;92(2):276-284.
6. Klaes R, Benner A, Friedrich T, et al. p16INK4a immunohistochemistry improves interobserver agreement in the diagnosis of cervical intraepithelial neoplasia. *Am J Surg Pathol.* 2002;26(11):1389-1399.
7. Negri G, Vitadello F, Romano F, et al. p16INK4a expression and progression risk of low-grade intraepithelial neoplasia of the cervix uteri. *Virchows Arch.* 2004;445(6):616-620.
8. Wentzensen N, von Knebel Doeberitz M. Biomarkers in cervical cancer screening. *Dis Markers.* 2007;23(4):315-330.
9. Cuschieri K, Wentzensen N. Human papillomavirus mRNA and p16 detection as biomarkers for the improved diagnosis of cervical neoplasia. *Cancer Epidemiol Biomarkers Prev.* 2008;17(10):2536-2545.
10. Tosunpol J, Arbyn M, Kyrgiou M, et al. p16INK4a immunostaining in cytological and histological specimens from the uterine cervix: a systematic review and meta-analysis. *Cancer Treat Rev.* 2009;35(3):210-220.
11. Dijkstra MG, Heideman DA, de Roy SC, et al. p16 (INK4a) immunostaining as an alternative to histology review for reliable grading of cervical intraepithelial lesions. *J Clin Pathol.* 2010;63(11):972-977.
12. Bergeron C, Ordi J, Schmidt D, et al. European CINtec Histology Study Group. Conjunctive p16INK4a testing significantly increases accuracy in diagnosing high-grade cervical intraepithelial neoplasia. *Am J Clin Pathol.* 2010;133:395-406.
13. Galgano MT, Castle PE, Atkins KA, et al. Using biomarkers as objective standards in the diagnosis of cervical biopsies. *Am J Surg Pathol.* 2010;34:1077-1087.
14. Massad LS, Einstein MH, Huh WK, et al. 2012 Updated consensus guidelines for the management of abnormal cervical cancer screening tests and cancer precursors. *J Low Genit Tract Dis.* 2013;17(5 Suppl 1):S1-S27.
15. Arbyn M, Kyrgiou M, Simoons C, et al. Perinatal mortality and other severe adverse pregnancy outcomes associated with treatment of cervical intraepithelial neoplasia: meta-analysis. *BMJ.* 2008;337:a1284.
16. Albrechtsen S, Rasmussen S, Thoresen S, et al. Pregnancy outcome in women before and after cervical conisation: population based cohort study. *BMJ.* 2008;337:a1343

17. Sadler L, Saftlas A, Wang W, et al. Treatment for cervical intraepithelial neoplasia and risk of preterm delivery. *JAMA*. 2004;291(17):2100-2106.
18. Park KJ, Soslow RA. Current concepts in cervical pathology. *Arch Pathol Lab Med*. 2009;133(5):729-738.
19. Stoler MH, Schiffman M. Atypical squamous cells of undetermined significance-low-grade squamous intraepithelial lesion triage study (ALTS) group. Interobserver reproducibility of cervical cytologic and histologic interpretations: realistic estimates from the ASCUS-LSIL Triage Study. *JAMA*. 2001;285(11):1500-1505.
20. Stoler MH, Ronnett BM, Joste NE, New Mexico HPV Pap registry steering committee, et al. The interpretive variability of cervical biopsies and its relationship to HPV status. *Am J Surg Pathol*. 2015;39(6):729-736.
21. Reuschenbach M, Wentzensen N, Dijkstra MG, et al. p16INK4a immunohistochemistry in cervical biopsy specimens: A systematic review and meta-analysis of the interobserver agreement. *Am J Clin Pathol*. 2014;142(6):767-72.
22. De Vet HCW, Knipschild PG, Schouten HJA, et al. Interobserver variation in histopathological grading of cervical dysplasia. *J Clin Epidemiol*. 1990;43:1395-8.
23. Creagh T, Bridger JE, Kupek E, et al. Pathologist variation in reporting cervical borderline epithelial abnormalities and cervical intraepithelial neoplasia. *J Clin Pathol*. 1995;48(1):p. 59-60.
24. Ceballos KM, Chapman DW, Daya D, et al. Reproducibility of the histological diagnosis of cervical dysplasia among pathologists from 4 continents. *Int J Gynecol Pathol*. 2008; 27(1):p. 101-7.
25. Guarch-Diaz CM, Suarez-Rincon AE, Vazquez-Camacho G, et al. p16INK4a immunohistochemistry improves the reproducibility of the histological diagnosis of cervical intraepithelial neoplasia in cone biopsies. *Gynecol Oncol*. 2008;111(1):120-124.
26. Horn LC, Reichert A, Oster A, et al. Immunostaining for p16INK4a used as a conjunctive tool improves interobserver agreement of the histologic diagnosis of cervical intraepithelial neoplasia. *Am J Surg Pathol*. 2008;32(4):502-512.
27. Sayed K, Korourian S, Ellison DA, et al. Diagnosing cervical biopsies in adolescents: the use of p16 immunohistochemistry to improve reliability and reproducibility. *J Low Genit Tract Dis*. 2007;11(3):141-146.
28. Darragh TM, Colgan TJ, Cox JT, et al. The lower anogenital squamous terminology standardization project for HPV-associated lesions: Background and consensus recommendations from the College of American Pathologists and the American Society for Colposcopy and Cervical Pathology. *J Low Genit Tract Dis*. 2012;16:205-242. Erratum in *J Low Genit Tract Dis*. 2013;17:36.
29. Stoler M, Bergeron C, Colgan TJ, et al. Tumors of the uterine cervix. In: Kurman RJ, Carcangiu ML, Herrington CS, et al. (Eds.), *WHO classification of tumors of female reproductive organs*. Lyon, France: IARC and WHO, 2014;169-2068.

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