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(...) = potenzialmente presente

Adattato da:


- 1 Stapleton, J.I. (1995). Host Immune Response to Hepatitis A Virus. *J Inf Dis* 171(suppl 1), 89-14.
- 2 Hollinger, F.B. et al. (2007). Hepatitis A virus. In: *Fields Virology: Knipe, D.M., Howley, P.M. (a cura di), 5ª ed.*, Lippincott Williams and Wilkins, Philadelphia, USA. Capitolo 27, 911-947.
- 3 Hadem, J. and Manns, M.P. (2007). Immune Response to Hepatitis A and E Viruses. Role in Disease Pathogenesis and Viral Elimination. In: Gershwin, M.E., Manns, M.P., Vierling, J.M., Springer, Link (Online service), editors. *Liver Immunology Principles and Practice*. Totowa, NJ: Humana Press Inc., 163-77.
- 4 Roque-Afonso, A.M. et al. (2010). Hepatitis A virus: serology and molecular diagnostics. *Future Virology* 5(2), 233-242.
- 5 Organizzazione Mondiale della Sanità (2011). The immunological basis for immunization series: module 18: hepatitis A. Disponibile all'indirizzo: https://apps.who.int/iris/bitstream/handle/10665/44570/9789241501422_eng.pdf?sequence=1. Accesso del 27 ottobre 2023
- 6 Salete de Paula, V. (2012). Laboratory diagnosis of hepatitis A. *Future Virol* 7(5), 461-472.
- 7 U.S. Department of Health and Human Services, Centers for Disease Control and Prevention (2015). Viral Hepatitis Serology Training, Hepatitis A. Disponibile all'indirizzo: <https://www.cdc.gov/hepatitis/resources/professionals/training/serology/training.htm>. Accesso del 27 ottobre 2023
- 8 Pischke, S. et al. (2018). Hepatitis A. In: Maus, S. et al. *Hepatology: A Clinical Textbook (9ª ed., pagg. 31-38)*. Medizin Fokus Verlag. Disponibile all'indirizzo: <https://www.hepatologytextbook.com/>. Accesso del 27 ottobre 2023

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
Adattato da:

- Hollinger, F.B. et al. (2007). Hepatitis A virus. In: *Fields Virology*. Knipe, D.M., Howley, P.M. (a cura di), 5a ed., Lippincott Williams and Wilkins, Philadelphia, USA. Capitolo 27, 911-947.
- Roque-Afonso, A.M. et al. (2010). Hepatitis A virus: serology and molecular diagnostics. *Future Virology* 5(2), 233-242
- Salete de Paula V (2012). Laboratory diagnosis of hepatitis A. *Future Virol* 7(5), 461-472.
- Ontario Association of Medical Laboratories (2010). Guidelines for Testing for Viral Hepatitis CLP12 and Interpretation of Viral Hepatitis Laboratory Test Results. Disponibile all'indirizzo: <https://oamtl.com/guidelines/>. Accesso del 27 ottobre 2023
- U.S. Department of Health and Human Services, Centers for Disease Control and Prevention (2015). Epidemiology and Prevention of Vaccine-Preventable Diseases - Hepatitis A. Disponibile all'indirizzo: <https://www.cdc.gov/vaccines/pubs/pinkbook/hepa.html>. Accesso del 27 ottobre 2023
- Public Health England. UK Standards for Microbiology Investigations - Hepatitis A virus acute infection serology (2019). *Virology* 27(4). Disponibile all'indirizzo: https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/772178/V_2714.pdf. Accesso del 27 ottobre 2023
- Oregon Health Authority. Public Health Division. Hepatitis A Investigative Guidelines (2019). Disponibile all'indirizzo: <https://www.oregon.gov/oha/PH/DISEASES/CONDITIONS/COMMUNICABLE/DISEASE/REPORTING/COMMUNICABLE/DISEASE/REPORTING/GUIDELINES/Documents/hepa.pdf>. Accesso del 27 ottobre 2023

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
² Elgouhari, H.M. et al. (2008). Hepatitis B virus infection: understanding its epidemiology, course, and diagnosis. *Cleve. Clin. J. Med.* 75, 881-889.

³ Fourati, S. and Pawlotsky, J.M. (2016). Recent advances in understanding and diagnosing hepatitis B virus infection. *F1000Res* 5:F1000 Faculty Rev-2243.

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


(...) = potenzialmente presente

Adattato da:

- 1 Elgouhari, H.M. et al. (2008). Hepatitis B virus infection: understanding its epidemiology, course, and diagnosis. *Cleve Clin J Med* 75, 881-889.
- 2 Fourati, S. and Pawlotsky, J.M. (2016). Recent advances in understanding and diagnosing hepatitis B virus infection. *F1000Res* 5:F1000 Faculty Rev-2243.
- 3 Lok, A.S. et al. (2017). Hepatitis B cure: From discovery to regulatory approval. *J Hepatol* 67, 847-861.
- 4 European Association for the Study of the Liver (EASL). 2017 Clinical Practice Guidelines on the management of hepatitis B virus infection. *J Hepatol* 67, 370-398.

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
Adattato da:

- 1 Elgouhari, H.M. et al. (2008). Hepatitis B virus infection: understanding its epidemiology, course, and diagnosis. *Cleve Clin J Med* 75, 881-889.
- 2 Fourati, S. and Pawlotsky, J.M. (2016). Recent advances in understanding and diagnosing hepatitis B virus infection. *F1000Res* 5:F1000 Faculty Rev-2243.
- 3 Lok, A.S. et al. (2017). Hepatitis B cure: From discovery to regulatory approval. *J Hepatol* 67, 847-861.
- 4 Davison, S.A. and Strasser S.I. (2014). Ordering and interpreting hepatitis B serology. *BMJ* 348, g2522.

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

Adattato da:

- Dufour, D.R. et al. (2001). Diagnosis and monitoring of hepatic injury. I. Performance characteristics of laboratory tests. *Clin Chem* 46(12), 2027-49.
- Hoofnagle, J.H. (2002). Course and Outcome of Hepatitis C. *Hepatology* 36, S21-S29.
- U.S. Department of Health and Human Services, Centers for Disease Control and Prevention (2015). Viral Hepatitis Serology Training, Hepatitis C. Disponibile all'indirizzo: <https://www.cdc.gov/hepatitis/resources/professionals/training/serology/training.htm>. Accesso del 27 ottobre 2023
- Linee guida dell'OMS sul test dell'epatite B e C. Geneva: World Health Organization; 2017 Feb. Fig. 1, Approximate time course for HCV virological and serological markers in chronic HCV infection. Disponibile all'indirizzo: <https://www.ncbi.nlm.nih.gov/books/NBK442283/figure/annex6.fig1/>. Accesso del 27 ottobre 2023
- Manns, M.P. et al. (2017). Hepatitis C virus infection. *Nat Rev Dis Prim* 3, 17006.
- Ahmad, J. (2017). Hepatitis C. *BMJ* 358, j2861.
- Austria, A. and Wu, G.Y. (2018). Occult Hepatitis C Virus Infection: A Review. *J Clin Transl Hepatol* 6, 155-160.

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
Elecsys® Immunoassays

Epatite virale	
Anti-HAV totale	
Anti-HAV IgM	
HBsAg	🔴
HBsAg	
HBsAg quantitativo	
Anti-HBs	🔴
Anti-HBc	🔴
Anti-HBc IgM	
Anti-HBe	
HBeAg quantitativo	
HBeAg	
Anti-HCV	🔴
Anti-HCV Duo	🔴

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
Adattato da:

- EASL Recommendations on Treatment of Hepatitis C (2018). *J Hepatol* 69, 461-511.
- Centers for Disease Control and Prevention (2013). Testing for HCV Infection: An Update of Guidance for Clinicians and Laboratorians. *MMWR* 62(18), 362-65.
- AASLD-ISA. Recommendations for testing, managing, and treating hepatitis C. <http://www.hcvguidelines.org>. Accessed 27/Oct/2023

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
Riferimenti:

- Aggarwal R, Goel A. Natural History, Clinical Manifestations, and Pathogenesis of Hepatitis E Virus Genotype 1 and 2 Infections. *Cold Spring Harbor Perspect Med.* 2019;9(7):a032136.
- Webb GW, Dalton HR. Hepatitis E: an underestimated emerging threat. *Ther Adv Infect Dis.* 2019;6:1-18.
- Lhomme S, et al. Screening, diagnosis and risks associated with Hepatitis E virus infection. *Exp Rev Anti-Inf Ther.* 2019;17:403-418.
- Kamar N, Izopet J, Pavio N, et al. Hepatitis E virus infection. *Not Rev Dis Primers.* 2017;3, 17086 (2017).
- OMS. <https://www.who.int/news-room/fact-sheets/detail/hepatitis-e>. Accesso del 5 febbraio 2024 Johns Hopkins Medicine. <https://www.hopkinsmedicine.org/health/conditions-and-diseases/hepatitis/hepatitis-a>. Accesso del 5 febbraio 2024
- Johns Hopkins Medicine. <https://www.hopkinsmedicine.org/health/conditions-and-diseases/hepatitis/hepatitis-a>. Accessed 5Feb2024

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
Riferimenti:

- 1 Kac P, Karra R. A Review of the Diagnosis and Management of Hepatitis E. *Curr Treat Options Infect Dis.* 2020;12:310-320.
- 2 Aslan AI, Balaban HY. Hepatitis E virus: Epidemiology, diagnosis, clinical manifestations, and treatment. *World J Gastroenterol.* 2020;26:5543-5560.
- 3 Public Health England. Public health operational guidelines for hepatitis E. Health protection response to reports of hepatitis E infection. 2019 Guidelines. Internet [aggiornato settembre 2019, citato maggio 2023]. Disponibile all'indirizzo: https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/845090/Public_Health_Operational_Guidelines_for_Hepatitis_E-protection_response.pdf.
- 4 Lhomme S, et al. Screening, diagnosis and risks associated with Hepatitis E virus infection. *Exp Rev Anti-Inf Ther.* 2019;17:403-418.
- 5 Associazione europea per lo studio del fegato (European Association for the Study of the Liver, EASL). Clinical Practice Guidelines on hepatitis E virus infection. *J Hepatol* 2018;68:1256-1271.

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
Adattato da:

- Aggarwal R, Goel A. Natural History, Clinical Manifestations, and Pathogenesis of Hepatitis E Virus Genotype 1 and 2 Infections. *Cold Spring Harbor Perspect Med.* 2019;9(7):a032136.
- Webb GW, Dalton HR. Hepatitis E: an underestimated emerging threat. *Ther Adv Infect Dis.* 2019;6:1-18.
- Lhomme S, et al. Screening, diagnosis and risks associated with Hepatitis E virus infection. *Exp Rev Anti-Infect Ther.* 2019;17:403-418.
- Kamar, N., Izopet, J., Pavio, N. et al. Hepatitis E virus infection. *Nat Rev Dis Primers.* 2017;3, 17086 (2017).
- OMS. <https://www.who.int/news-room/fact-sheets/detail/hepatitis-e>. Accesso del 5 febbraio 2024
- Johns Hopkins Medicine. <https://www.hopkinsmedicine.org/health/conditions-and-diseases/hepatitis/hepatitis-a>. Accesso del 5 febbraio 2024

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
Adattato da:

- 1 Abravanel F, et al. Diagnostic and management strategies for chronic hepatitis E infection. *Exp Rev Anti-Infect Ther*. 2023;21:143-148.
- 2 Kar, P, Karra R. A Review of the Diagnosis and Management of Hepatitis E. *Curr Treat Options Infect Dis*. 2020;12:310-320.
- 3 Aslan AI, Balaban HY. Hepatitis E virus: Epidemiology, diagnosis, clinical manifestations, and treatment. *World J Gastroenterol*. 2020;26:5543-5560.
- 4 Public Health England. Public health operational guidelines for hepatitis E. Health protection response to reports of hepatitis E infection. 2019 Guidelines. Internet [aggiornato settembre 2019, citato maggio 2023]. Disponibile all'indirizzo: https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/845090/Public_Health_Operational_Guidelines_for_Hepatitis_E-protection_response.pdf.
- 5 Lhomme S, et al. Screening, diagnosis and risks associated with Hepatitis E virus infection. *Exp Rev Anti-Infect Ther*. 2019;17:403-418.
- 6 Associazione europea per lo studio del fegato (European Association for the Study of the Liver, EASL). Clinical Practice Guidelines on hepatitis E virus infection. *J Hepatol*. 2018;68:1256-1271. Murali AR, et al. Chronic hepatitis E: A brief review. *World J Hepatol*. 2015;7:2194-201.
- 7 Murali AR, et al. Chronic hepatitis E: A brief review. *World J Hepatol*. 2015;7:2194-201.

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
Adattato da:

- 1 Prince, H.E. and Lapé-Nixon, M. (2014). Role of Cytomegalovirus (CMV) IgG Avidity Testing in Diagnosing Primary CMV Infection during Pregnancy. *Clin Vaccine Immunol* 21(10), 1377-1384.
- 2 Revello, M.G. and Gerna, G. (2002). Diagnosis and management of human cytomegalovirus infection in the mother, fetus, and newborn infant. *Clin Microbiol Rev* 15, 680-715.
- 3 Duff, P. (2010). Diagnosis and management of CMV Infection in Pregnancy. *Perinatology* 1, 1-6.
- 4 Centers for Disease Control (CDC) (2008). Knowledge and Practices of Obstetricians and Gynecologists Regarding Cytomegalovirus Infection During Pregnancy. *MMWR* 57(9), 65-68. Disponibile all'indirizzo: <https://www.cdc.gov/mmwr/preview/mmwrhtml/mm5703a2.htm> Accessed 27/01/2023
- 5 Davis, N.L. et al. (2017). Cytomegalovirus infection in pregnancy. *Birth Defects Research* 109, 336-346.

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Algoritmo di test¹⁻⁴


Adattato da:

- Munro, S.C. et al. (2005). Diagnosis of and screening for cytomegalovirus infection in pregnant women. *J Clin Microbiol* 43, 4713-4718.
- Duff, P. (2007). A thoughtful algorithm for the accurate diagnosis of primary CMV infection in pregnancy. *Am J Obstet Gynecol* 196, 196-197.
- Guerra, B. et al. (2007). Impact of diagnostic and confirmatory tests and prenatal counseling on the rate of pregnancy termination among women with positive cytomegalovirus immunoglobulin M antibody titers. *Am J Obstet Gynecol* 196, 221 e221-226.
- Lazzarotto, T. et al. (2004). Congenital cytomegalovirus infection: recent advances in the diagnosis of maternal infection. *Hum Immunol* 65, 410-415.

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2 Lambert, N. et al. (2015). Rubella. *Lancet* **385**, 2297-2307.

3 Vauloup-Fellous, C. and Grangeot-Keros, L. (2007). Humoral immune response after primary rubella virus infection and after vaccination. *Clin Vaccine Immunol* **14**, 644-647.

4 Abernathy, E. et al. (2009). Confirmation of rubella within 4 days of rash onset: comparison of rubella virus RNA detection in oral fluid with immunoglobulin M detection in serum or oral fluid. *J Clin Microbiol* **47**, 182-188.

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




Adattato da:

- 1 Riedel, A. et al. (2013). P5.071 Evaluation of Elecsys Immunoassay System for Determination of Type-Specific IgG Antibodies to HSV-1 and HSV-2. *Sex Transm Infect* 89(Suppl 1), A1-A428.
- 2 Patel, R. et al. (2001). European guideline for the management of genital herpes. *Int J STD AIDS* 12(Suppl 3), 34-39.
- 3 Morrow, R. and Friedrich, D. (2006). Performance of a novel test for IgM and IgG antibodies in subjects with culture-documented genital herpes simplex virus-1 or -2 infection. *Clin Microbiol Infect* 12, 463-469.
- 4 Whitley, R.J. and Miller, R.L. (2001). Immunologic approach to herpes simplex virus. *Viral Immunol* 14, 111-118.
- 5 Gardella, C. and Brown, Z.A. (2007). Managing genital herpes infections in pregnancy. *Cleve Clin J Med* 74, 217-224.
- 6 American Sexual Health Association. Herpes resource center: testing. Disponibile all'indirizzo: <http://www.ashsexualhealth.org/stdsstis/herpes/herpes-testing>. Accesso del 27 ottobre 2023
- 7 Sénat, M.V. et al. (2018). Prevention and management of genital herpes simplex infection during pregnancy and delivery: Guidelines from the French College of Gynaecologists and Obstetricians (FIGO). *Eur J Obstet Gynecol Repr Biol* 224, 93-101.
- 8 Workowski, K.A. et al. (2015). Sexually Transmitted Diseases Treatment Guidelines. *MMWR Recomm Rep* 64(3), 1-140.
- 9 Patel, R. et al. (2017). European guideline for the management of genital herpes. *Int J STD AIDS* 28(14), 1-14.
- 10 Groves, M.J. (2016). Genital Herpes: A Review. *Am Fam Physician* 93, 928-934.

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Algoritmo di test¹⁻⁴

Adattato da:

- Sénat, M.V. et al. (2018). Prevention and management of genital herpes simplex infection during pregnancy and delivery: Guidelines from the French College of Gynaecologists and Obstetricians (CNSOF). *Eur J Obstet Gynecol Repr Biol* 224, 93-101.
- Workowski, K.A. et al. (2015) Sexually Transmitted Diseases Treatment Guidelines, 2015. *MMWR Recomm Rep* 64(3), 1-140.
- Brown, Z.A. (2004). Use of Herpes Type-specific Serology to Prevent Neonatal Herpes Simplex Virus Infection. *NeoReviews* 5(1), e16-e21.
- Brown, Z.A. et al. (2005). Genital herpes complicating pregnancy. *Obstet Gynecol* 106, 845-856.

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

Adattato da:

1 Hess, R. (2004). Routine Epstein-Barr virus diagnostics from the laboratory perspective: still challenging after 35 years. *J Clin Microbiol* 42(8), 3381-7.

2 Middeldrop, J.A. (2015). Epstein-Barr virus-specific humoral immune responses in health and disease. In: C. Münz (ed.), *Epstein Barr Virus Volume 2, Current Topics in Microbiology and Immunology* 391, pagg. 289-322. Springer International Publishing Switzerland.

3 De Paschale, M. and Clerici, P. (2012). Serological diagnosis of Epstein-Barr virus infection: Problems and solutions. *World J Virol* 1(1), 31-43.

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
Adattato da:

- 1 De Paschale, M. and Cleirici, P. (2012). Serological diagnosis of Epstein-Barr virus infection: Problems and solutions. *World J Virol* 1(1), 31-43.
- 2 Public Health of England (PHE) (2019). UK Standards for Microbiology Investigations. Epstein-Barr virus serology. *Virology* 26(6), 2-8. Disponibile all'indirizzo: https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/773292/N_26i6.pdf. Ultimo accesso: [29 ottobre 2019].
- 3 Henle, G. et al. (1974). Antibodies to Epstein-Barr virus-associated nuclear antigen in infectious mononucleosis. *J Inf Dis* 130, 231-9.

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Epatite virale

Anti-HAV total	
Anti-HAV IgM	
HBsAg	●
HBsAg di conferma	
HBsAg quantitativo	
Anti-HBs	●
Anti-HBc	●
Anti-HBc IgM	
Anti-HBe	
HBeAg	
HBeAg quantitativo	
Anti-HCV	●
HCV Duo#	●

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© 2024 Roche
Roche Diagnostics S.p.A.