# Multicentre evaluation of a new strip-based blood glucose system (cobas<sup>®</sup> pulse, Roche Diagnostics) for near-patient testing in critical and non-critical care setting



Michael Goodman,<sup>1</sup> Nancy Breitenbeck,<sup>2</sup> William Clarke,<sup>3</sup> Cristina Canepa,<sup>5</sup> Edward Ganser,<sup>6</sup> Barbara Goldsmith,<sup>7</sup> Gary Headden,<sup>8</sup> Debra Hoppensteadt,<sup>9</sup> Birgit Klapperich,<sup>5</sup> Ryan Matika,<sup>10</sup> Gabrielle Miles,<sup>4</sup> Elizabeth Moser,<sup>4</sup> James H. Nichols,<sup>11</sup> Zubaid Rafique,<sup>12</sup> Ann Stankiewicz,<sup>4</sup> Jonathan Swanson,<sup>13</sup> Nam Tran,<sup>14</sup> Alan H. B. Wu,<sup>15</sup> Brad Karon<sup>16</sup>

<sup>1</sup>Department of Surgery, University of Cincinnati, CH, USA; <sup>2</sup>NB Research, Inc., Indianapolis, IN, USA; <sup>3</sup>Department of Pathology, Johns Hopkins University School of Medicine, Baltimore, MD, USA; <sup>4</sup>Roche Diagnostics Operations, Indianapolis, IN, USA; <sup>5</sup>Roche Diagnostics International AG, Rotkreuz, Switzerland; <sup>6</sup>Roche Diabetes Care Inc., Indianapolis, IN, USA; <sup>9</sup>Molecular Pharmacology and Neuroscience, Loyola University Chicago, IL, USA; <sup>10</sup>Department of Anesthesiology, and Immunology, Vanderbilt University Medical Center, Nashville, TN, USA; <sup>12</sup>Emergency Medicine, Baylor College of Medicine – Ben Taub Hospital, Houston, TX, USA; <sup>13</sup>Department of Pathology and Laboratory Medicine, University of California Davis Medical Center, Sacramento, CA, USA; <sup>15</sup>University of California, San Francisco, CA, USA; <sup>16</sup>Division of Clinical Core Laboratory Services, Department of Laboratory Medicine and Pathology, Mayo Clinic, Rochester, MN, USA.



### INTRODUCTION

- cobas<sup>®</sup> pulse (Roche Diagnostics GmbH, Mannheim, Germany) is a new near-patient whole blood (wb) glucose testing system for use in professional healthcare environments.
- No clinically relevant interference has been identified during rigorous evaluations of drugs, endogenous substances, naturally occurring sugars and artificial sweeteners.

 
 Table 2. Baseline characteristics of patients who contributed
arterial, venous and heel stick samples in the multicentre study

Charac-	Arterial	Venous	Heel stick	Total <sup>a</sup>
teristic	samples	samples	samples	(N=975)

Single-centre study: Capillary wb samples.

- cobas<sup>®</sup> pulse met both acceptance criteria (Table 5).
- NOVA StatStrip<sup>®</sup> did not meet the CLSI 12/12.5 criteria.

Table 5. Accuracy of capillary wb glucose measurements according to CLSI POCT12-A3 criteria

### **Breaknoint and criteria**

- The performance of this system has been evaluated in two separate studies that are the subject of this poster:
- A multicentre study that used arterial, venous and neonatal heel stick wb samples collected from patients in intensive critical care and noncritical health care settings.
- A single-centre study that used capillary, wb and finger stick samples collected from adult patients in a non-critical health care setting.

### METHODS

- The performance of the cobas<sup>®</sup> pulse system and Nova StatStrip<sup>®</sup> (Nova Biomedical, Waltham, MA, USA) were compared with a plasma-based hexokinase method on the Roche cobas<sup>®</sup> 6000 system (GLUC3).
- An aliquot from each sample (split sample) was measured on each of the three systems.
- All cobas<sup>®</sup> pulse measurements were performed by clinical staff who were self-trained on the use of the instrument using only the materials (package insert, quick start guide and operating manual) provided with the system.
- Accuracy was assessed by Clinical and Laboratory Standards Institute (CLSI) POCT12-A3 criteria (Table 1).
- Clinical significance was assessed by plotting pairs of results obtained with the cobas<sup>®</sup> pulse system and the hexokinase method on a Parkes error grid (Pfützner A. et al. J Diab Sci Technol 2013; 7: 1275-81).

 
 Table 1. CLSI POCT12-A3 criteria for point-of-care blood glucose
testing in acute and chronic care facilities

Breakpoint	100 mg/dL (5.55 mmol/L)	75 mg/dL (4.2 mmol/L)	
Threshold	12 mg/dL of the reference test result at <100 mg/dL, or 12.5% of the reference test result at ≥100 mg/dL	15 mg/dL of the reference test result at <75 mg/dL, or 20% of the reference test result at $\ge$ 75 mg/dL	
Acceptance	≥95% of results within threshold	≤2% of results exceed threshold	
Abbreviation used in this poster	12/12.5	15/20	

Gender				
Male	194 (61.4)	399 (56.5)	61 (50.8)	544 (55.8)
Female	122 (38.6)	307 (43.5)	59 (49.2)	431 (44.2)
Age				
≤28 days	n=62	n=16	n=120	198 (20.3)
29 days to <18 years	n=17	n=49	_	66 (6.8)
≥18 years	n=237	n=641	_	711 <sup>a</sup> (72.9)

### Values are n (%).

<sup>a</sup>Total unique patients (some patients contributed both arterial and venous samples).

In the single centre study, 702 capillary measurements were performed on samples collected from 117/130 enrolled non-critically ill patients at one US site (Figure 1B, Table 3).

 
 Table 3. Baseline characteristics of patients who contributed
capillary samples in the single-centre study

Characteristic	N (%)
Gender	
Male	45 (38.5)
Female	72 (61.5)
Age, years	
22 to 45	4 (3.5)
46 to 60	29 (24.8)
61 to 80	69 (59.0)
≥81	15 (12.8)
Diabetes	
Yes	99 (84.6)
No	18 (15.4)

	DIEARPUIIL AILU CILLEITA				
Within ! 0. (! 12 mg/c	(100 mg/dL) 67 mmol/L 1L)/12.5%, (%)	4.2 mmol/L (75 mg/dL) Outside by ! 0.83 mmol/L (! 15 mg/dL)/20%, n/N (%)			
cobas <sup>®</sup> pulse vs hexo-kinase	NOVA StatStrip <sup>®</sup> vs hexo-kinase	cobas <sup>®</sup> pulse vs hexo-kinase	NOVA StatStrip <sup>®</sup> vs hexo-kinase		
698/702 (99.4)	109/116 (94.0)	2/702 (0.3)	1/116 (0.9)		

Results in the shaded cell do not meet the acceptance criteria. n/N, number of samples meeting or not meeting acceptance criteria/ total number of samples.

Assessment of risk: Parkes error grid.

In the multicentre study (Figure 3), 99.4% (1135/1142) of the arterial, venous and heel stick wb glucose measurements performed on the cobas<sup>®</sup> pulse were within risk zone A, which corresponds to no effect on clinical action, whereas 0.6% of the results (7/1142) were in risk zone B, which represents altered clinical action of little or no effect on clinical outcome. No results were in zones C–E.

Figure 3. Parkes error grid of cobas<sup>®</sup> pulse system versus hexokinase test results for all arterial, venous and heel stick blood samples (N=1142) (multicentre study)



### RESULTS

### **Patients and samples.**

- In the multicentre study, 1290 wb samples were collected, of which 1142 were evaluable (316 arterial, 706 venous and 120 heel stick).
- These samples were collected from 975 patients (711 adult, 66 paediatric and 198 neonatal) ranging in age from <23 hours to >80 years, with and without diabetes in critical and non-critical care settings at 12 US sites (Figure 1A).

Multicentre study: Arterial, venous and neonatal heel stick wb samples. cobas<sup>®</sup> pulse fulfilled both acceptance criteria (Table 4).

- The cobas<sup>®</sup> pulse system also fulfilled both criteria in a sub-analysis that
- included arterial samples collected from neonates. NOVA StatStrip<sup>®</sup> did not meet the criteria for arterial, venous and heel stick
- samples at the 100 mg/dL breakpoint.
- NOVA StatStrip<sup>®</sup> also did not fulfil either criteria in the sub-analysis of neonatal samples.

Table 4. Accuracy of wb glucose measurements according to CLSI **POCT12-A3** criteria

	Breakpoint and criteria			
	5.55 mmol/L (100 mg/dL) Within ! 0.67 mmol/L (! 12 mg/dL)/12.5%, n/N (%)		4.2 mmol/L (75 mg/dL) Outside by ! 0.83 mmol/L (! 15 mg/dL)/20%, n/N (%)	
Sample type	cobas® pulse vs hexo- kinase	NOVA Stat-Strip® vs hexo- kinase	cobas® pulse vs hexo- kinase	NOVA Stat-Strip® vs hexo- kinase
Arterial (all)	312/316 (98.7)	277/316 (87.7)	2/316 (0.6)	3/316 (0.9)
Arterial (neo-nates)	61/62 (98.4)	45/62 (72.6)	1/62 (1.6)	2/62 (3.2)
Venous	691/706 (97.9)	652/705 (92.5)	8/706 (1.1)	6/705 (0.9)
Heel stick	114/120 (95.0)	111/120 (92.5)	0/120 (0)	2/120 (1.7)

In the single-centre study (Figure 4), 100% (702/702) of capillary blood glucose values were within risk zone A.

Figure 4. Parkes error grid of cobas<sup>®</sup> pulse system versus hexokinase test results for capillary blood samples (N=702) (single-centre study)





Characteristics of patients who contributed arterial, venous and heel stick samples are presented in Table 2.

Results in shaded cells do not meet the acceptance criteria. n/N, number of samples meeting or not meeting acceptance criteria/total number of samples.

### CONCLUSIONS

- The results of these two studies confirm the accuracy of the cobas<sup>®</sup> pulse system when used as intended in a diverse patient population by healthcare professionals working in hospital and clinic settings.
- The cobas<sup>®</sup> pulse system met CLSI POCT12-A3 accuracy criteria for arterial, venous, capillary and heel stick wb samples, and the performance exceeded that of NOVA StatStrip<sup>®</sup>.
- A clinical risk assessment of the data using Parkes error grids shows that the results obtained with the cobas<sup>®</sup> pulse system in these studies pose no clinically significant medical risks.

## DISCLOSURES

Dr. Michael Goodman reports participating in a Physician Advisory Board for Terumo BCT. All authors received research support in the form of medical writing assistance from Roche Diagnostics GmbH, Mannheim, Germany. Please see the abstract for all author disclosures. This study was funded by Roche Diagnostics.

## ACKNOWLEDGEMENTS

The authors would like to acknowledge the sustained collaborative effort and dedication of the entire Roche Point-of-Care Hospital Glucose Development Team in the realisation of this project, and Andre Schuetzenmeister, Manager Biostatistics, Biostatistics & Data Science, Roche Diagnostics GmbH, Penzberg, Germany for the Parkes error grid for the multicentre evaluation. Support for third-party medical writing assistance for this poster, furnished by Blair Jarvis, MSc, was provided by Roche Diagnostics GmbH, Mannheim, Germany.

> Copies of this poster obtained through quick response (QR) code are for personal use only and may not be reproduced without permission from the author of this poster. [DOI 10.1515/cclm-2021-5034]



This presentation is the intellectual property of the authors. Contact <u>goodmamd@ucmail.uc.edu</u> for permission to reprint, and/or distribute.

# **A New Digital Point-of-Care Tool With Advanced Blood Glucose Measuring Technology**

Simon Gessler, Dr. rer. nat.1 (D), Ulrich Porsch, MSc1 (D), David C. Klonoff, MD, FACP, FRCP (Edin), Fellow AIMBE2 (D), and Michael Marquant, Dipl.Ing. - Diplom Ingenieur der physikalischen Technik1

Journal of Diabetes Science and Technology 1–5 © 2022 Diabetes Technology Society

DOI: 10.1177/19322968221092763 journals.sagepub.com/home/dst

# ABSTRACT

cobas<sup>®</sup> pulse is a point-of-care blood glucose (BG) measuring system for multiple-patient use in professional healthcare settings. The system provides advances in connectivity and BG measuring technology, and has multiple fail-safes to improve accuracy and reduce the risk of user error. Flavin adenine dinucleotide-dependent glucose dehydrogenase on the working electrode catalyzes oxidation of β-D-glucose in the blood sample. A redox mediator/electron acceptor, on both the working and the counter electrode, facilitates diffusion of electrons in proportion to the glucose concentration and compensates for the effects of potential interfering agents. During development, >1 million test strip measurements were performed using >8000 test scenarios to refine the algorithm model. No clinically relevant interference was identified with extreme variations in blood properties and drugs in whole blood samples.



### **ELECTROCHEMISTRY OF GLUCOSE DETECTION IN THE SYSTEM**



Glucose detection occurs on gold electrodes in the reagent-coated slot dye. The working electrode is covered by a reagent layer (FADGDH) that catalyzes oxidation of  $\beta$ -D-glucose in the blood sample. On both, the working and the counter electrode, a redox mediator/electron acceptor facilitates diffusion of electrons from the co-factor to the electrode surface proportional to the concentration of glucose in the blood sample.

The FADGDH enzyme system is specific for glucose, and does not react with other saccharides in blood, which greatly reduces the risk for false positives.

### **COMPENSATION FOR ELECTRO-ACTIVE SUBSTANCES**

A key goal during the development of the system was absence of clinically relevant interferences. This was achieved by having the counter electrode generate a response current that closes the measurement circuit and by a change in polarity during the measurement sequence.

This process compensates for the effects of electro-active substances in blood samples.

\* PET: Polyethylenterephthalat; FADGH: flavin adenine dinucleotide-dependent glucose dehydrogenase.