



**From prolonged
time to diagnosis**



**to accelerated
diagnosis of AMI**

Faster, evidence-based algorithms
for treatment decisions based on the
2020 NSTEMI-ACS Guidelines from the
European Society of Cardiology

WHERE CARE LEADS



References

1. Collet, J.P. et al. (2021). Eur Heart J 42(14), 1289 – 367
2. Giannitsis, E. et al. (2021). Clin Res Cardiol, doi: 10.1007/s00392-021-01821-2
3. Ambavane, A. et al. (2017). PLoS One 12(11), e0187662
4. Thygesen, K. et al. (2018). Circulation 138(20), e618 – 51
5. Reichlin, T. et al. (2012) Arch Intern Med 172(16), 1211 – 8
6. Reichlin, T. et al. (2015) CMAJ 187(8), E243 – 52;
7. Mueller, C. et al. (2016) Ann Emerg Med 68(1), 76 – 87
8. Mokhtari, A. et al. (2016) JACC 67(13), 1531 – 40
9. Twerenbold, R. et al.(2019) J Am Coll Cardiol 74(4), 483 – 94
10. Chew, D.P. et al. (2019) Circulation 140(19), 1543 – 56
11. Shiozaki, M. et al. (2017) Int J Cardiol 249, 32 – 5
12. Stoyanov, K.M. et al. (2020). EHJ ACC 9(1), 39 – 51
13. Roffi, M. et al. (2016). Eur Heart J 37(3), 267 – 315
14. Elecsys Troponin T-hs and Troponin T hs STAT method sheets (2021).
Product numbers 09315322 190, 09315349 190, 09315357 190

Published by

Roche Diagnostics International Ltd
CH-6343 Rotkreuz
Switzerland

www.diagnostics.roche.com

© 2021

Biomarkers for the diagnosis of ACS



Key update from the
2020 NSTEMI-ACS
Guidelines from the ESC

1

Measurement of a **biomarker of cardiomyocyte** injury, preferably **hs-cTn**, is mandatory in all patients with suspected NSTEMI-ACS¹

2

If the clinical presentation is compatible with myocardial ischaemia, then a **dynamic elevation** of cTn >99th percentile of healthy individuals indicates Myocardial Infarction (MI)¹

3

The **ESC 0 /1-h algorithm** is now recommended as the first-choice diagnostic algorithm (Class I B)¹

4

As an alternative to the ESC 0/-1h algorithm, it is recommended to use the ESC **0/2-h algorithm**, if an hs-cTn test with a validated 0/2-h algorithm is available (Class I B)¹

5

The 0/3-h algorithm recommendation was **downgraded** to a Class IIa B recommendation¹


6

While cTnT-hs and hs-cTnI have comparable diagnostic accuracy, **cTnT-hs** has greater **prognostic** accuracy¹

7

Biomarkers should always be **used in combination** with **clinical assessment** and **12-lead ECG** in the diagnosis, risk stratification, and treatment of patients with suspected NSTEMI-ACS¹


The ESC 0 /1-h algorithm to accelerate the time to diagnosis and to drive immediate therapeutic consequences for the diagnosis of ACS



0/1-h algorithm benefits



Reduce length of observation time²



Reduced overall length of ED stay²




Lower hospital & overall AMI diagnostic cost^{2,3}

Risk of MI at index visit was <0.3% for the rule-out group, ~10% for the observe group, and >65% for the rule-in group; Risk of 30-day MACE was <0.5% for the rule-out group, 15–20% for the observe group, and >70% for the rule-in group.

* If chest pain onset was >3 h. AMI, acute myocardial infarction; CCTA, coronary computed tomography angiography; CCU, coronary care unit; ECG, electrocardiogram; ED, emergency department; ESC, European Society of Cardiology; MACE, major adverse cardiac event; MI, myocardial infarction; NSTEMI-ACS, non ST-elevation acute coronary syndrome

The first algorithm to rule-in or rule-out AMI within 0 to 1 h using cTnT-hs is confirmed by the 2020 ESC Guidelines


The ESC 0/1-h algorithm complemented with cut-offs specific for cTnT-hs¹



New in the guidelines: The ESC 0/2-h algorithm

Indicated as an alternative to the ESC 0/1-h algorithm

The ESC 0/2-h algorithm complemented with cut-offs specific for cTnT-hs¹



* Applicable for chest pain patients with onset longer than 3 h

AMI, acute myocardial infarction; ECG, electrocardiogram ESC, European Society of Cardiology

Caveats to consider with rapid algorithms¹



Algorithms should always be interpreted in conjunction with clinical and ECG findings¹




Troponin release is time-dependent. An additional cTn measurement at 3 h should be considered in patients presenting <1 h after chest pain onset and triaged towards rule-out¹



Late increases in cTn concentrations have been described. Serial cTn testing should be carried out if clinical suspicion remains high or in patients with recurrent chest pain¹

Assessing the rise and/or fall of cTn levels can help differentiate AMI from other conditions associated with cardiomyocyte injury⁴



¹Those with very low 0 h concentrations/chest pain onset >3 h
ESC, European Society of Cardiology; MI, myocardial infarction

Practical guidance for the implementation of the ESC 0/1 cTnT-hs algorithm¹



Collect blood samples for hs-cTn measurements at admission (0 h) and 1 h later¹



The team should not wait for other clinical details or pending results¹



The exact time of 0 h blood draw should be recorded to accurately estimate the time window (± 10 min) of the 1 h blood draw¹



In case the 1 h (± 10 min) collection is not possible, blood should be collected at 2 h and the 0/2-h algorithm should be used¹



Take note: Early measurements may seem unnecessary in some patients.* However, evaluating troponin dynamics early accelerates the process of ruling in and ruling out MI.

Timing of the blood draws and time of clinical decision differ

An important point to consider when using the 0/1-h algorithm



The ESC 2020 NSTEMI-ACS Guidelines recommend the measurement of cardiac troponins with high-sensitivity assays immediately after admission. The results should be available within 60 minutes of blood sampling (Class I, Level B)¹

The ESC 0/1-h algorithm is used irrespective of the local turnaround time. In the 0/1-h algorithm, 0 h and 1 h refer to the **timepoint at which blood is taken**. To calculate the earliest timepoint for clinical decision, the local TAT should be added to the time of blood collection¹



Compared to other algorithms, the 0/1-h algorithm possesses clinical and economic benefits, independent of the local TAT¹



Quicker diagnosis + shorter hospital stay = reduced resource utilisation³



Fewer blood draws, ECGs and imaging studies are required³

The performance and safety of rapid rule-out using the ESC 0/1-h algorithm with cTnT-hs is confirmed by multiple studies⁵⁻¹¹

Publications and trials using cTnT-hs values for patient assignment to the rule-out zone


Study	NPV (95% CI)	All-cause mortality or MACE in the rule-out zone
-------	--------------	--


0 h < 12 ng/L and Δ 1 h < 3 ng/L

APACE⁵	100%	30 days all-case mortality: 0.2% 2 years all-cause mortality: 1.9%
APACE⁶	99.9% (99.3–100%)	30 days all-case mortality: 0% 2 years all-cause mortality: 1.1%
TRAPID-AMI⁷	99.1% (98.2–99.7%)	30 days all-case mortality: 0.1% 2 years all-cause mortality: 0.7%
Mokhtari et al.⁸	NPV* for 30 days MACE: 99.5% (98.6–99.9%)*	30 days MACE*: 0.5% and 0% without UA


0 h < 5 ng/L or 0 h < 12 ng/L and Δ 1 h < 3 ng/L


APACE⁹	100%	30 days and 1-year all-case mortality: 0.2%
RAPID-TnT¹⁰	99.6% (99.0–99.9%) for 30 days death or MI	30 days all-case mortality and MI: 0.4%
Shiozaki et al.¹¹	100% (96.8–100%)	30 days all-case mortality: 0%

 **cTnT-hs can be used as an aid for early discharge and out-patient management for patients suspected of ACS**

 The **high NPV** (99.1–100%) of the 0/1-h algorithm and the **low 30-day mortality** (0–0.4%) in the rule-out zone **confirms** the **safety** of this approach for **early discharge**

In three prospective studies and more than 8000 patients interventional studies confirm the excellent performance of the 0/1-h algorithm using cTnT-hs^{9,10,12}




 The ESC 0/1-h cTnT-hs algorithm triages over three-quarters of ED chest pain patients: **75–84% of patients** triaged within 1 hour (plus lab time)

The safety of the ESC 0/1-h algorithm using cTnT-hs is confirmed by the low incidence of 30-day MACE or mortality in rule-out patients*

Twerenbold R et al.⁹
30-day MACE rate:
0.2% in the rule-out group
0.1% among patients triaged to outpatient care[†]

Chew DP et al.¹⁰
30-day all-cause death or MI rate:
0.5% in the rule-out group
0.3% among patients discharged directly from the ED

Stoyanov KM et al.¹²
30-day mortality rate:
0.4% in the rule-out group
0.08% among patients triaged to outpatient care[†]


 The **low 30-day MACE** rate (0.1–0.3%) among patients discharged directly from the ED confirms the safety of this approach for **early discharge**

* When using the extended algorithm adding ECG and patient history as recommended by the ESC guidelines. ACS, acute coronary syndrome; CI, confidence interval; ECG, electrocardiogram; NPV, negative predictive value; MACE, major adverse cardiac event; MI, myocardial infarction; UA, unstable angina

* Rule-out group indicates patients recommended to be ruled-out for MI by the algorithm; † Outpatient care indicates patients in which the final management decision made was direct discharge from the ED CI, confidence interval; ED, emergency department; ESC, European Society of Cardiology; MACE, major adverse cardiac events; MI, myocardial infarction

The performance of the ESC 0/1-h algorithm using cTnT-hs is robust in important patient subgroups (e.g. early presenters, sex, age)⁹

Forest plots illustrating efficacy and safety of the ESC 0/1-h cTnT-hs algorithm in pre-defined subgroups⁹




30-day MACE rates were low in all subgroups

For early presenters (CPO ≤3 h), the 0/1-h algorithm triaged **80%** to rule-out or rule-in, with a **0%** 30-day MACE rate in rule-out patients⁹


The ESC 0/1-h algorithm using cTnT-hs reduces LoS in ED & acute care^{9,10,12}

Median reduction of 1–2 hours of length of stay* compared with standard of care^{9,10,12}



Improves ED patient flow

- Triages **75–84%** of patients to rule-in OR rule-out within 1 hour (plus lab time)^{9,10,12}
- Rule-out of AMI possible in **62–72%**^{9,10,12} of patients with a NPV of **99.6%**⁹



Shortens length of stay (compared to 0/3-h algorithm) without compromising safety

- Significantly **increases ED discharge** rates compared with the standard of care^{9,10,12}
- Reduces ED LoS by 1.0 – 2.25 h^{9,10,12}
- Reduces acute care LoS by 1.1 h¹²

*Comparisons are with pre implementation / standard algorithm. AMI, acute myocardial infarction; ED, emergency department; ESC, European Society of Cardiology; LoS, length of stay; NPV, negative predictive value



Conclusions from the three studies^{9,10,12}

The implementation of the ESC 0/1-h algorithm using cTnT-hs¹³




cTnT-hs is safe

- Low 30-day event rate 0.2 - 0.5% (MACE, all cause mortality and MI, or mortality)^{9,10,12}
- 0.08% 30-day mortality in patients triaged to out-patient care¹²



cTnT-hs is feasible

- 94% algorithm adherence⁹
- 45 minute median reduction in the time between initial and second cTnT-hs sample¹²



cTnT-hs does not increase the use of diagnostic resources^{9,10,12}

- Significant reduction in functional cardiac testing¹⁰
- No increase in coronary angiography^{10,12}

Biomarkers as tools for prognostication and risk assessment




Beyond its diagnostic role, it is recommended to measure hs-cTn serially for the estimation of prognosis **(Class I B)**¹



cTnT-hs and hs-cTnI possess comparable diagnostic accuracy, however, **cTnT-hs has greater prognostic accuracy**¹



hs-cTn **adds prognostic information** with regards to short- and long-term mortality to clinical and ECG variables¹



The higher the hs-cTn levels, the greater the risk of death. Serial measurements can be used to **detect peak levels** of cardiac troponin and risk stratify patients with established MI¹



BNP/NT-proBNP, serum creatinine and eGFR also affect prognosis and should be **measured concomitantly**¹

Conclusions




The 0/1-h algorithm (best option, blood draw at 0h and 1h) and the 0/2-h algorithm (second best option, blood draw at 0h and 2 h) are the preferred rapid algorithms and are recommended with a **Class IB** recommendation¹



Both algorithms have **data driven assay-specific cut-offs**, derived and validated to achieve predefined performance characteristics, for both rule-out and rule-in of AMI, justifying immediate therapeutic consequences¹



Multiple hs-cTn assays can be used with both algorithms and are supported by real life studies documenting safe and efficient implementation¹



The 0/1-h algorithm using cTnT-hs is **validated in numerous prospective trials**^{9,10,12}



Roche cTnT-hs is the first hs-cTn assay to introduce a **rapid rule-out claim for early discharge** and **out-patient management** for patients suspected of ACS¹⁴

Elecsys® Troponin T hs intended use¹⁴

Supported by numerous (> 1300) publications and has multiple intended uses



1. Aid in diagnosis

Aid in the differential diagnosis of acute coronary syndrome to identify necrosis, e.g. acute myocardial infarction¹⁴



5. Selection of therapy and intervention

The test may also be useful for the selection of more intensive therapy and intervention in patients with elevated levels of cardiac troponin T¹⁴




2. Rapid rule-out

Aid for early discharge and out-patient management for patients suspected of acute coronary syndrome (ACS)¹⁴



6. Perioperative use

To predict the perioperative risk of major adverse cardiac events and for the aid in diagnosis of (MINS) and PMI in non-cardiac surgeries¹⁴




3. Risk stratification - ACS

Risk stratification of patients presenting with acute coronary syndrome¹⁴



7. General population

To aid in stratifying the long-term risk of cardiovascular death, myocardial infarction, coronary revascularization, heart failure, or ischemic stroke and all-cause mortality in asymptomatic individuals¹⁴



4. Risk stratification – chronic renal failure

Cardiac risk stratification in patients with chronic renal failure¹⁴