

# cobas p 480 Instrument

## *Automating your primary vial preprocessing steps*

The **cobas p 480** Instrument improves laboratory efficiency by allowing valuable technician time to be used more productively, eliminating repetitive, manual sample handling, improving workflow and reducing risk of contamination, human error and workplace injury.



### Improving laboratory efficiency

- Accepts PreservCyt®, SurePath™ liquid based cytology vials as well as **cobas**® PCR Media and **cobas**® PCR Cell Collection Media primary vials
- Processes four vials simultaneously
- Intuitive interface requires minimal training
- High throughput automation allows a single **cobas p 480** Instrument to support up to two **cobas**® 4800 Systems, delivering



Up to 376 **cobas**® 4800 System results in 8.7 hours



Up to 470 **cobas**® 4800 System results in an 8.5 shift\*

\*Final run completes overnight

Decapping workflow using 1 **cobas p 480** Instrument and 2 **cobas**® 4800 Systems

### Reduces hands on time and repetitive motions with four unique workflows




#### Decapping

Removes caps from primary tubes for testing on the **cobas**® 4800 System



#### Recapping

Recaps sample vials with new caps to avoid contamination



#### Aliquoting

Aliquots configurable volumes from PreservCyt®, SurePath™ and **cobas**® PCR Cell Collection Media primary vials into barcode matched secondary tubes

Compatible with SurePath™ vials with plastic inserts



#### Reagent addition and heating\*

Addition of **cobas**® Sample Prep Buffer and incubation of specimens collected in SurePath™ medium to reverse cross-linking and free nucleic acids for testing<sup>1</sup>



\***cobas**® Sample Prep Buffer and use of SurePath™ Medium with the **cobas**® HPV Test are not approved by the FDA for use in the United States.



*Life needs answers*

### Improves sample reproducibility and process reliability

- Sample chain of custody is assured with primary and secondary vial barcode matching
- All vials are spun prior to opening to remove potentially contaminating droplets from sample caps
- Precision pipetting using CO-RE tip, Total Aspirate and Dispense Monitoring and Anti-Droplet Control technologies reduce opportunities for contamination and ensure sample integrity
- No LIS or data connection required
- Printable reports capture all sample ID's, sample error and reagent lot and expiration information



### System Specifications

<b>Power</b>	
Power consumption	Max. 1000 VA
Line Voltage	115 VAC (-15%) to 230 (+10%) VAC
Line frequency	50/60 Hz $\pm$ 5
Power interruption	115 VAC: 10 A 230 VAC: 5 A
<b>Interfaces</b>	
Instrument to PC	USB
<b>Dimensions</b>	
HxWxD (with full cover)	90.5cm x 166.5cm x 101cm (35.6" x 65.55" x 39.76")
<b>Weight</b>	~180 kg

<sup>1</sup>Kiernan, J. A. Preservation and retrieval of antigens for immunohistochemistry – methods and mechanisms Part 2. Retrieving masked antigens  
Department of Anatomy and Cell Biology. The Cutting Edge: 5-11. <http://publish.uwo.ca/~jkiernan/FixAnti2.pdf>

COBAS, COBAS P and LIFE NEEDS ANSWERS are trademarks of Roche.  
All other product names and trademarks are the property of their prospective owners.

© 2015 Roche Molecular Systems, Inc.

Roche Molecular Diagnostics  
4300 Hacienda Drive  
Pleasanton, CA 94588  
USA  
[www.molecular.roche.com](http://www.molecular.roche.com)

