 FOUNDATIONONE® CDx

 FOUNDATIONONE® LIQUID CDx

 FOUNDATIONONE® HEME

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
FOUNDATIONONE® CDx

PATIENT: XXXX TUMOR TYPE: Ovary high grade serous carcinoma COUNTRY CODE: CH REPORT DATE: XXXX

ABOUT THE TEST: FoundationOne® CDx is a next-generation sequencing (NGS) based assay that identifies genomic findings within hundreds of cancer-related genes.

PATIENT: DISEASE: Ovary high grade serous carcinoma NAME: XXXX DATE OF BIRTH: XXXX1955 SEX: Female MEDICAL RECORD #: XXXX

PHYSICIAN: ORDERING PHYSICIAN: USZ, Pathologie MEDICAL FACILITY: Institut fuer Pathologie und Molekularpathologie ADDITIONAL RECIPIENT: None MEDICAL FACILITY ID: 201897 PATHOLOGIST: .

SPECIMEN: SPECIMEN SITE: Spleen SPECIMEN ID: XXXX SPECIMEN TYPE: FPE DATE OF COLLECTION: Invalid date SPECIMEN RECEIVED: XXXX

Report Highlights

- Targeted therapies with NCCN categories of evidence in this tumor type: **Niraparib** (p. 8), **Olaparib** (p. 9), **Rucaparib** (p. 10)
- Evidence-matched clinical trial options based on this patient's genomic findings: (p. 12)

Genomic Signatures

Loss of Heterozygosity score - 32.0%
 Microsatellite status - MS-Stable
 Tumor Mutational Burden - 1 Muts/Mb

Gene Alterations

For a complete list of the genes assayed, please refer to the Appendix.

BARD1 deletion exon 3
 TP53 R158G
 BRAF amplification
 KEL amplification

2 Disease relevant genes with no reportable alterations: **BRCA1, BRCA2**

1 See About the Test in appendix for details.

GENOMIC SIGNATURES	THERAPIES WITH CLINICAL RELEVANCE (IN PATIENT'S TUMOR TYPE)	THERAPIES WITH CLINICAL RELEVANCE (IN OTHER TUMOR TYPE)
Loss of Heterozygosity score = 32.0%	Niraparib [ZA]	Talazoparib [6]
10 Trials see p. 12	Olaparib [7]	
Microsatellite status = MS-Stable	Rucaparib [10]	
Tumor Mutational Burden = 1 Muts/Mb	No therapies or clinical trials. see Genomic Signatures section	
GENE ALTERATIONS	THERAPIES WITH CLINICAL RELEVANCE (IN PATIENT'S TUMOR TYPE)	THERAPIES WITH CLINICAL RELEVANCE (IN OTHER TUMOR TYPE)
BARD1 - deletion exon 3	none [8]	none
10 Trials see p. 14		
TP53 - R158G	none [11]	none
2 Trials see p. 16		

© 2022 Foundation Medicine, Inc. All rights reserved. PAGE 1 OF 23

Lab Analysis and Sequencing: University Hospital Zurich, Schmelzbergstrasse 12, 8091 Zurich, Switzerland
 Computational Analysis and Report Curation: 150 Second St, 1st Floor, Cambridge, MA 02141 / CLIA: 22D2027531

2

3

1

4

5

10

7

6

8

11

FOUNDATIONONE® CDx

PATIENT: XXXX TUMOR TYPE: Ovary high grade serous carcinoma COUNTRY CODE: CH REPORT DATE: XXXX

USZ# XXXX

TABLE THERAPEUTIC OR CLINICAL TRIAL OPTIONS

Biological and clinical significance, including prognostic, diagnostic, germline, and potential chemosensitivity data, are provided in the following section.

p. 7 **KEL - amplification** p. 7

With activity of certain drugs approved by applicable regulatory authorities (for example, the FDA, EMA, or country specific regulatory authorities) however, the agency for the patient's tumor type. Neither the therapeutic agents nor the trials identified are ranked in order of potential or predicted efficacy for this patient, nor are they ranked for the therapy section. The report should only be regarded and used as a supplementing source of information. All treatment decisions remain the full and various reasons further explained below both the therapies and the clinical trials listed in this report they not be complete and exhaustive. All drugs displayed in the only in the patient's tumor type. Please find the entire Suite prescribing information to www.f1m5med.com.

© 2022 Foundation Medicine, Inc. All rights reserved. PAGE 2 OF 23

Lab Analysis and Sequencing: University Hospital Zurich, Schmelzbergstrasse 12, 8091 Zurich, Switzerland
 Computational Analysis and Report Curation: 150 Second St, 1st Floor, Cambridge, MA 02141 / CLIA: 22D2027531

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Références:

1. FoundationOne®Heme Technical Specifications, 2021. www.foundationmedicine.com/genomic-testing/foundation-one-heme (dernière consultation: février 2022). 2. FoundationOne®CDx Technical Specifications, 2021. www.foundationmedicine.com/test/foundationone-cdx (dernière consultation: février 2022). 3. FoundationOne®LiquidCDx Technical Specifications, 2021. <https://www.foundationmedicine.com/test/foundationone-liquid-cdx> (dernière consultation: février 2022). 4. Frampton, GM et al. "Development and validation of a clinical cancer genomic profiling test based on massively parallel DNA sequencing." Nature biotechnology vol. 31,11 (2013): 1023-31. 5. Drilon, A et al. "Broad, Hybrid Capture-Based Next-Generation Sequencing Identifies Actionable Genomic Alterations in Lung Adenocarcinomas Otherwise Negative for Such Alterations by Other Genomic Testing Approaches." Clinical cancer research : an official journal of the American Association for Cancer Research vol. 21,16 (2015): 3631-9. 6. He, J et al. "Integrated genomic DNA/RNA profiling of hematologic malignancies in the clinical setting." Blood vol. 127,24 (2016): 3004-14. 7. FoundationOne®CDx: Technical Information Version 02. <https://www.rochefoundationmedicine.com/flcdxtech> (dernière consultation: février 2020). 8. Site internet de l'Hôpital universitaire de Zurich Pathologie moléculaire www.usz.ch/fachbereich/pathologie-molekularpathologie/angebot/molekulares-tumorprofilung (dernière consultation: février 2022). 9. Ross, Jeffrey S et al. "Comprehensive Genomic Profiling of Carcinoma of Unknown Primary Site: New Routes to Targeted Therapies." JAMA oncology vol. 1,1 (2015): 40-49. 10. Approbation de FoundationOne®Liquid CDx par la FDA, 2020. Disponible sous: <https://www.foundationmedicine.com/press-releases/445c1f9e-6cbb-488b-84ad-5f133612b721> (dernière consultation: mars 2022). 11. Planchard, D et al. "Correction to: "Metastatic non-small cell lung cancer: ESMO Clinical Practice Guidelines for diagnosis, treatment and follow-up"." Annals of oncology: official journal of the European Society for Medical Oncology vol. 30,5 (2019): 863-870. 12. NCCN Clinical Practice Guidelines in Oncology, Non-Small Cell Lung Cancer, Version 6.2020, June 2020 Disponible sous: https://www.nccn.org/professionals/physician_gls/pdf/nscl.pdf (dernière consultation: août 2020). 13. Cardoso, F et al. "Early breast cancer: ESMO Clinical Practice Guidelines for diagnosis, treatment and follow-up." Annals of oncology: official journal of the European Society for Medical Oncology vol. 30,10 (2019): 1674. 14. NCCN Clinical Practice Guidelines in Oncology, Breast Cancer, Version 5.2020, July 2020. Disponible sous: https://www.nccn.org/professionals/physician_gls/pdf/breast.pdf (dernière consultation: août 2020). 15. NCCN Clinical Practice Guidelines in Oncology, Prostate Cancer, Version 2.2021, May 2020. Disponible sous: https://www.nccn.org/professionals/physician_gls/pdf/prostate.pdf (dernière consultation: août 2021). 16. NCCN Clinical Practice Guidelines in Oncology, Ovarian Cancer, Version 1.2020, March 2020. https://www.nccn.org/professionals/physician_gls/pdf/ovarian.pdf (dernière consultation: août 2020). 17. Colombo, N et al. "ESMO-ESGO consensus conference recommendations on ovarian cancer: pathology and molecular biology, early and advanced stages, borderline tumours and recurrent disease." Annals of oncology: official journal of the European Society for Medical Oncology vol. 30,5 (2019): 672-705. 18. Parker, C et al. "Prostate cancer: ESMO Clinical Practice Guidelines for diagnosis, treatment and follow-up." Annals of oncology: official journal of the European Society for Medical Oncology vol. 31,9 (2020): 1119-1134. 19. Leibel, S et al. "Waiting for a breast biopsy. Psychosocial consequences and coping strategies." Journal of psychosomatic research vol. 55,5 (2003): 437-43. 20. Hayes Balmadrid, MA et al. "Anxiety prior to breast biopsy: Relationships with length of time from breast biopsy recommendation to biopsy procedure and psychosocial factors." Journal of health psychology vol. 22,5 (2017): 561-571. 21. Francis G, Stein S. "Circulating Cell-Free Tumour DNA in the Management of Cancer." International journal of molecular sciences vol. 16,614122-42. 22. Siravegna, G et al. "Clonal evolution and resistance to EGFR blockade in the blood of colorectal cancer patients." Nature medicine vol. 21,7 (2015): 795-801. 23. Remon, J et al. "Osimertinib benefit in EGFR-mutant NSCLC patients with T790M-mutation detected by circulating tumour DNA." Annals of oncology: official journal of the European Society for Medical Oncology vol. 28,4 (2017): 784-790. 24. Bidard, Francois-Clement et al. "Going with the flow: from circulating tumor cells to DNA." Science translational medicine vol. 5,207 (2013): 207ps14. 25. Villafior, Victoria et al. "Biopsy-free circulating tumor DNA assay identifies actionable mutations in lung cancer." Oncotarget vol. 7,41 (2016): 66880-66891. 26. Allen, Justin M et al. "Genomic Profiling of Circulating Tumor DNA in Relapsed EGFR-mutated Lung Adenocarcinoma Reveals an Acquired FGFR3-TACC3 Fusion." Clinical lung cancer vol. 18,3 (2017): e219-e222. 27. Scherer, F. "Capturing tumor heterogeneity and clonal evolution by circulating tumor DNA profiling." Tumor Liquid Biopsies (2020): 213-230. 28. Tsui, DWY et al. "Tumor fraction-guided cell-free DNA profiling in metastatic solid tumor patients." Genome medicine vol. 13,196. 31 May, 2021. 29. Taylor, BS et al. "Advances in sarcoma genomics and new therapeutic targets." Nature reviews. Cancer vol. 11,8 541-57. 14 Jul. 2011. 30. Cancer Genome Atlas Research Network. Cancer Genome Atlas Research Network. "Comprehensive and Integrated Genomic Characterization of Adult Soft Tissue Sarcomas." Cell vol. 171,4 (2017): 950-965.e28. 31. Boddü, Spandana et al. "Clinical Utility of Genomic Profiling in the Treatment of Advanced Sarcomas: A Single-Center Experience." JCO precision oncology vol. 2 (2018): 1-8. 32. Dong, L et al. "Clinical Next Generation Sequencing for Precision Medicine in Cancer." Current genomics vol. 16,4 (2015): 253-63. 33. NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines). https://www.nccn.org/professionals/physician_gls/default.aspx (dernière consultation: mars 2022).

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