

cobas[®] infinity laboratory solution 3.0

Product information

cobas[®] infinity laboratory solution

LIS/Host Communication

cobas[®] infinity supports standard communication protocols to connect with Laboratory Information Systems and other host systems.

Supported messages include:

- Order Query
- Demographic Updates
- Patient Results
- QC Results

- HL7 v2.x & v3.x
- Native XML
- SOAP Web Services
- IHE standard frameworks for Laboratory Testing Workflow (LTW)
- Laboratory Analytical Workflow (LAW)

Instrument Connectivity

cobas[®] infinity supports standard communication protocols to connect with laboratory instruments.

Supported messages include:

- Order Query
- Demographic Updates
- Patient Results
- QC Results

- Pre-analytics:
 - cobas[®] connection modules (CCM)
 - cobas[®] 8100 automated workflow series
 - MODULAR *PRE-ANALYTICS* (MPA)
 - cobas p 512, cobas p 612, cobas p 471, cobas p 671, cobas p 312 (incl. Dynamic Interface)
 - Other 3rd party systems
- Post-analytics:
 - cobas p 501/701
 - Rotastock
- Connectivity with 3rd party applications via (ASTM, HL7 v2/v3)

Server operating systems

(administrator rights required)

- Microsoft Windows Server 2008 R2 SP1 Enterprise Edition
- Microsoft Windows Server 2012 R2 Standard Edition
- Microsoft Windows Server 2016 Standard Edition

Server environment

Database server: InterSystem[®] Healthshare Foundation[®] Database.

Web Server: Apache[™] server.

Any virtualized system with resources dedicated following the outlined hardware specifications can be used as a virtualized environment.

Antivirus usage is allowed.

High availability solution options

Overview of the Central Lab deployment options (mirroring is the option recommended by InterSystems):

- **Mirroring (synchronous and asynchronous)**

Database mirroring is a high-availability solution that is used to increase database availability. In the event of a disaster, a standby copy of the database can be brought online.

- **Clustering**

Configuration based on a system configured in cluster. It consists of a set of disks (disk rack), shared by two identical servers, with a system of automatic switching between servers; it also has a RAID 5 for automatic load balancing on the disks.

Server hardware requirements

Size	S	M	L	XL	
Tubes/day (including aliquotes)	5,000	10,000	20,800	52,400	
Results/day	35,100	71,550	224,000	560,000	
Peak Tubes/hour (including reflex/rerun)	500	1,000	2,500	6,600	
Peak results/hour	3,510	7,155	26,923	70,534	
Users' Concurrency	8	15	30	50	
Hardware recommendations					
RAM (DDR4 - GB)	8–16	8–16	16–32	16–32	
CPU (cores)	4–6	6–8	20–24	24–32	
Storage (retention days)					
OS (GB)	36	36	36	36	
DB (RAID 10) SSD** (GB)	Data Traces	(0.5 GB x retention days) x 2 (0.5 GB x log days) x 2	(1 GB x retention days) x 2 (1 GB x log days) x 2	(3.5 GB x retention days) x 2 (2.7 GB x log days) x 2	(9 GB x retention days) x 2 (7 GB x log days) x 2
Images (RAID 10) SSD**	See estimating additional disk space requirements section				
Journals (RAID 10) SSD** (GB)	(1 GB x journal days) x 2	(15 GB x journal days) x 2	(50 GB x journal days) x 2	(90 GB x journal days) x 2	
Backups (RAID 5) HDD (2 days)* (GB)	DB x 1.2 x backup days	DB x 1.2 x backup days x 2	DB x 1.2 x backup days	DB x 1.2 x backup days	
Reports (RAID 5) HDD*	See estimating additional disk space requirements section				

* No fast disk needed ** SSD disk recommended. This value has been multiplied per 2 due to the RAID 10

Languages and regions

cobas® infinity laboratory solution supports multiple languages and regions (international time zones)

User management

- Support LDAP (Lightweight Directory Access Protocol)
- Support of level 3 of LOPD (Law of Personal Data Protection)

Security

- The use of firewalls and SSL encoding is recommended
- Protection against direct browsing is enabled

User interface

100% browser based user interface on the workstations (clients)

- Supported browsers*
 - Microsoft Edge 42
 - Google Chrome 69
 - Microsoft Internet Explorer® 11
 - Mozilla® Firefox® 60 (ESR version) and 38 (ESR version)
 - Client screen resolution – 1280 x 1024
 - Adobe® Flex® plug-in may be required for applications or special graphic software
- * Supported browsers will be updated with each version of cobas® infinity Laboratory Solutions

cobas® infinity IT solutions client PC requirements

Processor	RAM	Hard Disk	Network	CD-ROM	UPS	Monitor
Pentium Core 2 Duo 2.4 GHz	4 GB	1 HDD 80 GB	Megabit Ethernet	DVD drive	An uninterruptible power supply installation is recommended	At least 1280x1024 pixels with 65,536 colors (min.)

Printers and barcode printers (peripherals)

- Flexible customization of reports (DataVision, FOP) and scheduling of print jobs
- Print compatibility with any printer, optical pen, barcode reader and document scanner
- The following barcode printers are supported up-to-date:
 - Eltron TL27xx
 - Easy Coder 91 and Easy Coder C4
 - Intermec 3500, PC43D and compatible models
 - Zebra SP 10-z, 2746e, TLP 2844-Z, GX 420t and compatible models

Remote service and support

- SIP Service Intelligent Platform (non customer-facing tool that continuously monitors the health and versions of your customer software solution)
- Axeda

Additional drivers may be available upon request.
For a complete list of drivers, contact your local Roche representative.

Find out more on
diagnostics.roche.com