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LITERATURA

1. Syria Laperche, Nadine Le Marrec, Annie Girault, Francoise Bouchardeau, Annabelle Servant-Delmas, Michèle Maniez-Montreuil, Pierre Gallian, Thierry Levayer, Pascal Morel, and Nicole Simon. Simultaneous Detection of Hepatitis C Virus (HCV) Core Antigen and Anti-HCV Antibodies Improves the Early Detection of HCV Infection; *JOURNAL OF CLINICAL MICROBIOLOGY*, Aug. 2005, p. 3877-3883, Vol. 43, No. 8.
2. J. Morgan Freiman, MD, Trang M. Tran, Samuel G Schumacher, MSc, PhD, Laura F. White, PhD, Stefano Ongareello, PhD, Jennifer Cohn, MD, MPH, Philippa J. Easterbrook, MD, MPH, Benjamin P. Linas, MD, MPH, and Claudia M. Denking, MD, PhD. HCV Core Antigen Testing for Diagnosis of HCV Infection: A systematic review and meta-analysis; *Ann Intern Med.* 2016 September 06; 165(5): 345-355. doi:10.7326/M16-0065.
3. Cresswell, Martin Fisher, Daniel J. Hughes, Simon G. Shaw, Gary Homer, and Mohammed O. Hassan-Ibrahim. Hepatitis C Core Antigen Testing: A Reliable, Quick, and Potentially Cost-effective Alternative to Hepatitis C Polymerase Chain Reaction in Diagnosing Acute Hepatitis C Virus Infection; *Clinical Infectious Diseases** 2015; 60(2): 263-6.
4. Shah O. Dinesh, Chi D. Chang, Lily X. Jiang, Kevin Y. Cheng, A. Scott Muerhoff, Robin A. Gutierrez, Thomas P. Leary, Suresh M. Desai, Irene V. Batac-Herman, Vince A. Salbilla, Alla S. Haller, James L. Stewart, and George J. Dawson. Combination HCV core antigen and antibody assay on a fully automated chemiluminescence analyzer *Transfusion* 2003, 43 (8): 1067-74
5. Tuke PW, Grant PR, Waite J, Kitchen AD, Eglin RP, Tedder RS (2008). Hepatitis C virus window-phase infections: closing the window on hepatitis C virus. *Transfusion* 48: 594-600. Volume 43, August 2003 *TRANSFUSION* 1067, From the New Assay Development and Infectious Diseases; R & D, Abbott Diagnostics Division, Abbott Laboratories, Abbott Park, Illinois.
6. Roth WK, et al. International survey on NAT testing of blood donations: expanding implementation and yield from 1999 to 2009. *Vox Sang* 2012; 102: 82-90 [PMID: 21933190 DOI: 10.1111/j.1423-0410.2011.01506.x
7. Ross RS, Viazov S, Salloum S, Hilgard P, Gerken G, Roggendorf M. Analytical performance characteristics and clinical utility of a novel assay for total hepatitis C virus core antigen quantification. *J Clin Microbiol* 2010; 48: 1161-1168 [PMID: 20107102 DOI: 10.1128/JCM.01640-09].
8. Hans L Tillmann Hepatitis C virus core antigen testing: Role in diagnosis, disease monitoring and treatment; *World J Gastroenterol* 2014 June 14; 20(22): 6701-6706 ISSN 1007-9327 (print) ISSN 2219-2840 (online) *World J Gastroenterol* 2014 March 21; 20(11): 2948-2954 ISSN 1007-9327.
9. Neelam Marwaha, Suchet Sachdev. Current testing strategies for hepatitis C virus infection in blood donors and the way forward; *WJG 20th Anniversary Special Issues (2): Hepatitis C virus.*
10. Micha C. Nübling, Gabriele Unger, Michael Chudy, Steven Raia, Johannes Löwer Sensitivity of HCV core antigen and HCV RNA detection in the early infection phase; *Transfusion* vol. 42; August 2002, 1037-45.
11. Sebastiaan J. Hullege, Corine H. Geurtsvan-Kessel, Annemiek A. van der Eijk, Chris Ramakers and Bart J. A. Rijnders. HCV antigen instead of RNA testing to diagnose acute HCV in patients treated in the Dutch Acute HCV in HIV Study; *Journal of the International AIDS Society* 2017, 20:21621. <http://www.jiasociety.org/index.php/jias/article/view/21621> | <http://dx.doi.org/10.7448/IAS.20.1.21621>
12. Cornelia Ottiger*, Nicole Gygli, Andreas R. Huber. Detection limit of architect hepatitis C core antigen assay in correlation with HCV RNA, and renewed confirmation algorithm for reactive anti-HCV samples; *Journal of Clinical Virology* 58 (2013) 535-540.
13. Harun Khan, Andrew Hill, Janice Main, Ashley Brown, and Graham Cooke. Can Hepatitis C Virus Antigen Testing Replace Ribonucleic Acid Polymearse Chain Reaction Analysis for Detecting Hepatitis C Virus? A Systematic Review; *Open Forum Infectious Diseases**© The Author 2017. DOI: 10.1093/ofid/ofw252.
14. AK Reddy, KV Dakshinamurty, V Lakshmi. Utility of HCV core antigen elisa in the screening for hepatitis C virus infection in patients on hemodialysis *Indian J Med Microbiol.* 2006 Jan;24(1):c 55-7.
15. Peterson J, Green G, Iida K, et al. Detection of hepatitis C core antigen in the antibody negative "window" phase of hepatitis C infection. *Vox Sang* 2000; 78: 80-5.
16. Shaheen A.M, Amin A.M, Tawfik A.K, Ezz Elden M.M. Analytical evaluation of Hepatitis C Virus Core Antigen Detection as a Diagnostic Test in Diagnosis of HCV Infection; *Egyptian Journal of Medical Microbiology*, October 2014 Vol. 23, No. 41.
17. Katja Semea, Mario Poljaka, Dunja Z. Babic, Tina Mocilnik, Adriana Vince. The role of core antigen detection in management of hepatitis C: a critical review. *Journal of Clinical Virology* 32 (2005), 92-101.
18. Steven H. Kleinman, Nico Lelie, and Michael P. Busch. REVIEW – Infectivity of human immunodeficiency virus-1, hepatitis C virus, and hepatitis B virus and risk of transmission by transfusion. *TRANSFUSION* Volume 49, November 2009, 2454-2489.
19. 2017 European Guideline for the screening, prevention and initial management of hepatitis B & C infections in sexual health settings (2017 IUSTI Europe Hepatitis Guideline), <https://www.iusti.org/regions/europe/pdf/2017/IUSTIEuropeHepatitis2017.pdf>
20. Yuksel P, Caliskan R, Ergin S, Aslan M, Celik DG, Saribas S, Ziver T, Yalciner A, Kocazeybek B. New approaches to in vitro diagnosis of hepatitis C infection a reason for post transfusion hepatitis: Diagnostic value of determination of hepatitis C virus core antigen. *Transfus Apher Sci.* 2011 Dec;45(3): 247-50. doi: 10.1016/j.transci.2011.10.002. Epub 2011 Nov 2.
21. Salih Cetiner, Alev Cetin Duran, Filiz Kibar, Akgün Yaman Performance comparison of new generation HCV core antigen test versus HCV RNA test in management of hepatitis C virus infection. *Transfusion and Apheresis Science* 56 (2017) 362-366.
22. Kuo Yuan-Hung, Kuo-Chin Chang, Jing-Houng Wang, Pei-Shan Tsai, Shu-Feng Hung, Chao-Hung Hung, Chien-Hung Chen, and Sheng-Nan Lu. Is Hepatitis C Virus Core Antigen an Adequate Marker for Community Screening? *Journal of Clinical Microbiology* June 2012 Volume 50 Number 6, p. 1989-1993.
23. Turek P., Masopust. J. Činnost transfuzní služby v České republice v období 1989–2013. *Transfuzie Hematol. dnes*, 20, 2014, No. 4, p. 125-135.
24. APASL consensus statements and recommendation on treatment of hepatitis C - *Hepatology* Int. DOI 10.1007/s12072-016-9717-6.
25. AASLD: Testing, Evaluation, and Monitoring of Hepatitis C. From www.HCVGuidance.org on November 27, 2018.
26. European Association for the Study of the Liver. EASL Recommendations on Treatment of Hepatitis C 2018. *J Hepatol* (2018), <https://doi.org/10.1016/j.jhep.2018.03.026>
27. https://talk.ictvonline.org/ictv_wikis/flaviviridae/w/sg_flavi/56/hcv-classification