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## LITERATURA

1. Verlohren S, Stepan H, Dechend R (2012) Angiogenic growth factors in the diagnosis and prediction of pre-eclampsia. *Clin Sci* 122: 43–52.
2. Abalos E, Cuesta C, Carroli G, Qureshi Z, Widmer M, Vogel JB, Souza JP. Preeclampsia, eclampsia and adverse maternal and perinatal outcomes: a secondary analysis of the World Health Organization Multicountry Survey on Maternal and Newborn Health. *BJOG* 2014; 121 Suppl 1: 14–24. DOI: 10.1111/1471-0528.12629.
3. Duley L. The global impact of pre-eclampsia and eclampsia. *Semin Perinatol* 2009; 33: 130–137. DOI: 10.1053/j.semperi.2009.02.010.
4. American College of Obstetricians and Gynecologists, Task Force on Hypertension in Pregnancy. Hypertension in pregnancy (2013) Report of the American College of Obstetricians and Gynecologists' Task Force on hypertension in pregnancy. *Obstet Gynecol* 122(5): 1122–1131. <https://doi.org/10.1097/01.AOG.0000437382.03963.88>
5. Poon LC, Sheman A, Hyett JA et al. The International Federation of Gynecology and Obstetrics (FIGO) initiative on pre-eclampsia: a pragmatic guide for first-trimester screening and prevention. *Int J Gynaecol Obstet.* 2019; 145: 1–33.
6. Brown MA, Magee LA, Kenny LC, Karumanchi SA, McCarthy FP, Saito S et al (2018) Hypertensive disorders of pregnancy: ISSHP classification, diagnosis, and management recommendations for international practice. *Hypertension* 72: 24–43.
7. Brunelli VB, Prefumo F. Quality of first trimester risk prediction models for pre-eclampsia: a systematic review. *BJOG* 2015; 122: 904–914. DOI: 10.1111/1471-0528.13334.
8. Poon LC, Kametas NA, Maiz N, Akolekar R, Nicolaides KH. First-trimester prediction of hypertensive disorders in pregnancy. *Hypertension* 2009; 53: 812–818. DOI:10.1161/hypertensionaha.108.127977.
9. Tan MY, Syngelaki A, Poon LC, Rolnik DL, O'Gorman N, Delgado JL, Akolekar R, Konstantinidou L, Tsavdaridou M, Galeva S, Ajdacka U, Molina FS, Persico N, Jani JC, Plasencia W, Greco E, Papaioannou G, Wright A, Wright D, Nicolaides KH. Screening for pre-eclampsia by maternal factors and biomarkers at 11–13 weeks' gestation. *Ultrasound Obstet Gynecol* 2018; 52: 186–195. DOI: 10.1002/uog.19112.
10. Rolnik DL, Wright D, Poon LCY, Syngelaki A, O'Gorman N, de Paco Matallana C, Akolekar R, Cicero S, Janga D, Singh M, Molina FS, Persico N, Jani JC, Plasencia W, Papaioannou G, Tenenbaum-Gavish K, Nicolaides KH. ASPRE trial: performance of screening for preterm pre-eclampsia. *Ultrasound Obstet Gynecol* 2017; 50: 492–495. DOI: 10.1002/uog.18816.
11. Rolnik DL, Wright D, Poon LC, O'Gorman N, Syngelaki A, de Paco Matallana C, Akolekar R, Cicero S, Janga D, Singh M, Molina FS, Persico N, Jani JC, Plasencia W, Papaioannou G, Tenenbaum-Gavish K, Meiri H, Gizurarson S, Maclagan K, Nicolaides KH. Aspirin versus placebo in pregnancies at high risk for preterm preeclampsia. *N Engl J Med* 2017; 377: 613–622. DOI: 10.1056/NEJMoa1704559.
12. Wright D, Tan MY, O'Gorman N, Poon LC, Syngelaki A, Wright A, Nicolaides KH. Predictive performance of the competing risk model in screening for preeclampsia. *Am J Obstet Gynecol* 2019; 220: 199.e191–199.e113. DOI: 10.1016/j.ajog.2018.11.1087.
13. Zeisler H, Llorba E, Chantraine F, Vatish M, Staff AC, Sennstrom M, Olovsson M, Brennecke SP, Stepan H, Allegranza D, Dinkel C, Schoedl M, Dilba P, Hund M, Verlohren S. Soluble fms-Like Tyrosine Kinase-1-to-Placental Growth Factor Ratio and Time to Delivery in Women With Suspected Preeclampsia. *Obstet Gynecol* 2016; 128: 261–269. DOI: 10.1097/aog.0000000000001525.
14. Poon, L. C., Galindo, A. , Surbek, D. , Chantraine, F. , Stepan, H. , Hyett, J. , Tan, K. H. and Verlohren, S. (2019), From first-trimester screening to risk stratification of evolving pre-eclampsia in the second and third trimesters of pregnancy: a comprehensive approach. *Ultrasound Obstet Gynecol.* Accepted Author Manuscript. doi:10.1002/uog.21869.
15. Maynard SE, Min J-Y, Merchan J, Lim K-H, Li J, Mondal S et al (2003) Excess placental soluble fms-like tyrosine kinase 1 (sFlt1) may contribute to endothelial dysfunction, hypertension, and proteinuria in preeclampsia. *J Clin Invest* 111:649–658.
16. Levine RJ, Maynard SE, Qian C, Lim K-H, England LJ, Yu KF et al (2004) Circulating angiogenic factors and the risk of preeclampsia. *N Engl J Med* 350: 672–683.
17. Levine RJ, Lam C, Qian C, Yu KF, Maynard SE, Sachs BP et al (2006) Soluble endoglin and other circulating antiangiogenic factors in preeclampsia. *N Engl J Med* 355: 992–1005.
18. Verlohren S, Galindo A, Schlembach D, Zeisler H, Herraiz I, Moertl MG et al (2010) An automated method for the determination of the sFlt-1/PIGF ratio in the assessment of preeclampsia. *Am J Obstet Gynecol* 202: 161.e1–161.e11.
19. Verlohren S, Herraiz I, Lapaire O, Schlembach D, Moertl M, Zeisler H et al (2012) The sFlt-1/PIGF ratio in different types of hypertensive pregnancy disorders and its prognostic potential in preeclamptic patients. *Am J Obstet Gynecol* 206: 58.e1–58.e8.
20. Verlohren S, Herraiz I, Lapaire O, Schlembach D, Zeisler H, Calda P et al (2014) New gestational phase-specific cutoff values for the use of the soluble fms-like tyrosine kinase-1/placental growth factor ratio as a diagnostic test for preeclampsia. *Hypertension* 63: 346–352.
21. Zeisler H, Llorba E, Chantraine F, Vatish M, Staff AC, Sennström M et al (2016) Soluble fms-like tyrosine kinase-1-to-placental growth factor ratio and time to delivery in women with suspected preeclampsia. *Obstet Gynecol* 128: 261–269.
22. Zeisler H, Llorba E, Chantraine FJ, Vatish M, Staff AC, Sennström M et al (2018) The sFlt-1/PIGF Ratio: ruling out pre-eclampsia for up to 4 weeks and the value of retesting. *Ultrasound Obstet Gynecol* 53(3): 367–375. <https://doi.org/10.1002/uog.19178>.