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5. Schindewolf E, Khalek N, Johnson MB, et al. Expanding the fetal phenotype: Prenatal sonographic findings and perinatal outcomes in a cohort of patients with a confirmed 22q11.2 deletion syndrome. *Am J Med Genet Part A*. 2018; 176(8): 1735-1741. doi:10.1002/ajmg.a.38665.
6. Morrow BE, McDonald-McGinn DM, Emanuel BS, Vermeesch JR, Scambler PJ. Molecular genetics of 22q11.2 deletion syndrome. *Am J Med Genet Part A*. 2018; 176(10): 2070-2081. doi:10.1002/ajmg.a.40504.
7. Rauch A, Hoyer J, Guth S, et al. Diagnostic yield of various genetic approaches in patients with unexplained developmental delay or mental retardation. *Am J Med Genet Part A*. 2006; 140(19): 2063-2074. doi:10.1002/ajmg.a.31416.
8. Cheung ENM, George SR, Andrade DM, Chow EWC, Silversides CK, Bassett AS. Neonatal hypocalcemia, neonatal seizures, and intellectual disability in 22q11.2 deletion syndrome. *Genet Med*. 2014; 16(1): 40-44. doi:10.1038/gim.2013.71.
9. Bassett AS, McDonald-McGinn DM, Devriendt K, et al. Practical guidelines for managing patients with 22q11.2 deletion syndrome. *J Pediatr*. 2011; 159(2). doi:10.1016/j.jpeds.2011.02.039
10. Campbell IM, Sheppard SE, Crowley TB, et al. What is new with 22q? An update from the 22q and You Center at the Children's Hospital of Philadelphia. *Am J Med Genet Part A*. 2018; 176(10): 2058-2069. doi:10.1002/ajmg.a.40637.
11. Bevilacqua, E. et al. (2020) Performance of a targeted cell-free DNA prenatal test for 22q11.2 deletions in a large clinical cohort. Poster presented at the International Society of Ultrasound in Obstetrics and Gynecology's Virtual World Congress.
12. Horowitz A, Shifman S, Rivlin N, Pisanté A, Darvasi A (2005). „A survey of the 22q11 microdeletion in a large cohort of schizophrenia patients“. *Schizophr. Res.* 73 (2-3): 263-7. doi:10.1016/j.schres.2004.02.008.
13. Kagan K, Hoopmann M, Pfaff T et al., First Trimester Screening for Common Trisomies and Microdeletion 22q11.2 Syndrome Using Cell-Free DNA: A Prospective Clinical Study, *Fetal Diagn Ther* 2020; 47: 841-851.
14. Demografická ročenka České republiky, D.21, ČSÚ, roky 1989-2019.
15. Centre NC. Multiple pregnancy. *Heal (San Fr)*. 2011; September).
16. ACOG. Practice Bulletin 169: Multifetal gestations: Twin, triplet and higher-order multifetal pregnancies. *Obstet Gynecol*. 2016; 128(4): e131-e146.
17. Oepkes D, Suetters M. Antenatal fetal surveillance in multiple pregnancies. *Best Pract Res Clin Obstet Gynaecol*. 2017; 38: 59-70.
18. Khalil A, Rodgers M, Baschat A, et al. ISUOG Practice Guidelines: role of ultrasound in twin pregnancy. *Ultrasound Obstet Gynecol*. 2016; 47(2): 247-263.
19. American College of Obstetricians and Gynecologists, Society for Maternal-Fetal Medicine. *Obstetric Care Consensus No. 6*. *Obstet Gynecol*. 2017; 130(4): e187-e199.
20. Audibert F, Gagnon A. No. 262-Prenatal Screening for and Diagnosis of Aneuploidy in Twin Pregnancies. *J Obstet Gynaecol Canada*. 2017; 39(9): e347-e361.
21. Gil M, Accurti V, Santacruz B, Plana M, Nicolaides K. Analysis of Cell-Free DNA in Maternal Blood in Screening For Aneuploidies: Updated Meta-Analysis. *Ultrasound Obstet Gynecol*. 2017; 50: 302-314.
22. Struble C, Syngelaki A, Oliphant A, Song K, Nicolaides KH. Fetal Fraction Estimate in Twin Pregnancies Using Directed Cell-Free DNA Analysis. *Fetal Diagn Ther*. December 2013.
23. Canick JA, Kloza EM, Lambert-Messerlian GM, et al. DNA sequencing of maternal plasma to identify Down syndrome and other trisomies in multiple gestations. *Prenat Diagn*. 2012; 32(8): 730-734.
24. Grömminger S, Yagmur E, Erkan S, et al. Fetal Aneuploidy Detection by Cell-Free DNA Sequencing for Multiple Pregnancies and Quality Issues with Vanishing Twins. *J Clin Med*. 2014; 3(3): 679-692.
24. del Mar Gil M, Quezada MS, Bregant B, Syngelaki A, Nicolaides KH. Cell-Free DNA Analysis for Trisomy Risk Assessment in First-Trimester Twin Pregnancies. *Fetal Diagn Ther*. 2014; 35(3): 204-211.
26. Gil M, Galeva S, Jani J, Konstantinidou L, Akolekar R, Plana M, Nicolaides K. Screening for trisomies by cfDNA testing of maternal blood in twin pregnancy: update of The Fetal Medicine Foundation results and meta-analysis *Ultrasound Obstet Gynecol*. 2019 Jun; 53(6): 734-742.
27. Jones KJ, Wang E, Bogard P, et al. Performance of targeted cell-free DNA (cfDNA) analysis with microarray quantitation for assessment of fetal sex and sex chromosome aneuploidy risk. *Ultrasound Obstet Gynecol*. November 2017. doi:10.1002/uog.18968.
28. Stokowski R, Wang E, White K, et al. Clinical performance of non-invasive prenatal testing (NIPT) using targeted cell-free DNA analysis in maternal plasma with microarrays or next generation sequencing (NGS) is consistent across multiple controlled clinical studies. *Prenat Diagn*. 2015; 35: 1243-1246.
29. Bevilacqua E, Gil MM, Nicolaides KH, et al. Performance of screening for aneuploidies by cell-free DNA analysis of maternal blood in twin pregnancies. *Ultrasound Obstet Gynecol*. 2015; 45(1): 61-66.
30. Sarno L, Revello R, Hanson E, Akolekar R, Nicolaides KH. Prospective screening for trisomies by cell-free DNA testing of maternal blood in first trimester twin pregnancies. *Ultrasound Obstet Gynecol*. March 2016.
31. Galeva S, Konstantinidou L, Gil MM, Akolekar R, Nicolaides KH. Routine first-trimester screening for fetal trisomies in twin pregnancy: cell-free DNA test contingent on results from combined test. *Ultrasound Obstet Gynecol*. 2019 Feb; 53(2): 208-213.
32. Galeva S, Gil MM, Konstantinidou L, Akolekar R, Nicolaides KH. First-trimester screening for trisomies by cfDNA testing of maternal blood in singleton and twin pregnancies: factors affecting test failure. *Ultrasound Obstet Gynecol*. 2019 Jun; 53(6): 804-809.
33. Benn P, Borrell A, Chiu RWK, et al. Position statement from the Chromosome Abnormality Screening Committee on behalf of the Board of the International Society for Prenatal Diagnosis. *Prenat Diagn*. 2015; 35(8): 725-734.
34. Schmid M, White K, Stokowski R, Miller D, Bogard PE, Valmeekam V, Wang E. Accuracy and reproducibility of fetal-fraction measurement using relative quantitation at polymorphic loci with microarray. *Ultrasound Obstet Gynecol*. 2018 Jun; 51(6): 813-817.
35. Sparks AB, Struble CA, Wang ET, Song K, Oliphant A. Noninvasive prenatal detection and selective analysis of cell-free DNA obtained from maternal blood: evaluation for trisomy 21 and trisomy 18. *Am J Obstet Gynecol*. 2012; 206(4): 319.e1-9.
36. Brar H, Wang E, Struble C, Musci TJ, Norton ME. The fetal fraction of cell-free DNA in maternal plasma is not affected by a priori risk of fetal trisomy. *J Matern Fetal Neonatal Med*. 2013; 26(2): 143-145.
37. Wang E, Batey A, Struble C, Musci T, Song K, Oliphant A. Gestational age and maternal weight effects on fetal cell-free DNA in maternal plasma. *Prenat Diagn*. 2013; 33(7): 662-666.
38. Nicolaides KH, Musci, TJ, Struble, C Syngelaki, A Gil, M M. Assessment of Fetal Sex Chromosome Aneuploidy Using Directed Cell-Free DNA Analysis *Fetal Diagn Ther*. 2014; 35(1): 1-6.30.
39. Revello R, Sarno L, Ispas A, Akolekar R, Nicolaides KH. Screening for trisomies by cell-free DNA testing of maternal blood: consequences of failed result. *Ultrasound Obstet Gynecol*. January 2016.
40. ACOG. Screening for fetal aneuploidy. *Practice Bulletin Number 163*. *Obstet Gynecol*. 2016.
41. Bevilacqua E, Chen K, Wang Y, Doshi J, White K, de Marchin J, Conotte S, Jani JC, Schmid M. Cell-free DNA analysis after reduction in multifetal pregnancies. *Ultrasound Obstet Gynecol*. 2019 Jun 10.
42. Thurik FF, Ait Soussan A, Bossers B, et al. Analysis of false-positive results of fetal RHD typing in a national screening program reveals vanishing twins as potential cause for discrepancy. *Prenat Diagn*. April 2015.
43. Kelley JF, Henning G, Ambrose A, Adelman A. Vanished Twins and Misdiagnosed Sex: A Case Report with Implications in Prenatal Counseling Using Noninvasive Cell-Free DNA Screening. *J Am Board Fam Med*. 29(3): 411-413.
44. Vlková B, Hodosy J. Vanishing twin as a potential source of bias in non-invasive fetal sex determination: A case report. *J Obstet Gynaecol Res*. 2014; 40(4): 1128-1131.
45. Judah H, Gil MM, Syngelaki A et al. Cell-free DNA testing of maternal blood in screening for trisomies in twin pregnancy: cohort study at 10-14 weeks and updated meta-analysis, *Ultrasound Obstet Gynecol* 2021.
46. Norton E, Jacobsson B, Swamy G, et al., Cell-free DNA analysis for noninvasive examination of trisomy, *N Engl J Med*. 2015 Apr 23; 372(17): 1589-97.