



cobas c 701 and c 702 analyzers
Within-Run Precision Guidelines



COBAS and LIFE NEEDS ANSWERS are trademarks of Roche.

©2012-2020 Roche

Roche Diagnostics
9115 Hague Rd
Indianapolis, IN 46256
USA
www.roche.com
www.dialog.roche.com

The contents of this document, including all graphics and photographs, are the property of Roche. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Roche.

Every effort has been made to ensure that the information contained in this manual is accurate at the time of printing. Not all functionality described in this manual may be available to all users. Roche Diagnostics reserves the right to make any further required changes to software without prior notice. Such changes may not immediately be reflected in this document.

Any screenshots in this publication have been added exclusively for the purpose of illustration. Configurable and variable data such as parameters, results, path names etc. visible therein must not be used for laboratory purposes.

This document is intended for the US market only.

Caution: Federal law restricts this device to sale by or on the order of a physician

This document was created by the Roche Diagnostics Engineering Operations department. Direct questions or concerns regarding the contents of this document to:

Roche Diagnostics Corporation
Engineering Operations Department
9115 Hague Road
Indianapolis, IN 46256
USA

COBAS, COBAS C and LIFE NEEDS ANSWERS are trademarks of Roche.
All other trademarks are the property of their respective owners.

©2010-2020, Roche Diagnostics. All right reserved.

[VV-06833-02](#)

Distribution in USA by: Roche Diagnostics, Indianapolis, IN

This document is available on the Roche Diagnostics USA website at www.dialog.roche.com

Revision History

Revisions to this document are provided by Roche when necessary. No part of this document may be reproduced in any form or by any means without prior written consent.

Publication Reference Number	Date	Revision purpose
4157-00-1110 Version 0	November 2010	Document creation with cobas c 701 analyzer data
4157-01-1014 Version 1	October 2014	Added new cobas c 701 and c 702 data.
4157-02-0817 Version 2	August 2017	Added new HDLC4 and LDLC3 within run precision guidelines.
4157-03-0817 Version 3	August 2017	Updated Ion Selective Electrodes section.
VV-06833-02 Version 4	November 2020	Updated customer website URL; Deleted HDL-C plus gen.3, LDL-C plus gen.2, T4 and T Uptake; Added Ethyl Alcohol, B2-Microglobulin (urine), Tina-quant C-Reactive Protein IV and Transferrin ver2 (urine); Updated Ammonia gen.2, CK, Total Protein U/CSF (CSF), Total Protein U/CSF (urine), B2-Microglobulin (serum/plasma), Rheumatoid Factors II and Transferrin ver2 (serum/plasma).

How to use these guidelines:

Where the SD and CV% information is provided, the SD specification is used to determine pass/fail criteria if the actual sample concentration is below the stated sample concentration in the table. If the actual sample concentration is equal to or above the stated sample concentration in the table, the CV% specification is used to determine pass/fail criteria.

Assay	SD	Sample Concentration	%CV Guideline
Albumin BCG, gen. 2	≤ 0.07 g/dL	3.5 g/dL	≤ 2 %
Albumin BCP	≤ 0.07 g/dL	3.5 g/dL	≤ 2 %
Alk Phos/IFCC gen. 2	≤ 2 U/L	100 U/L	≤ 2 %
ALT	≤ 0.9 U/L	30 U/L	≤ 3 %
ALT/P-5-P	≤ 1.1 U/L	35 U/L	≤ 3 %
Ammonia gen. 2	≤ 4.3 µg/dL	85 µg/dL	≤ 5 %
p-Amylase/EPS (serum)	≤ 1 U/L	50 U/L	≤ 2 %
p-Amylase/EPS (urine)	≤ 11 U/L	350 U/L	≤ 3 %
t-Amylase/EPS (serum)	≤ 2 U/L	100 U/L	≤ 2 %
t-Amylase/EPS STAT (serum)	≤ 2 U/L	100 U/L	≤ 2 %
t-Amylase/EPS (urine)	≤ 14 U/L	460 U/L	≤ 3 %
t-Amylase/EPS STAT (urine)	≤ 14 U/L	460 U/L	≤ 3 %
AST	≤ 0.9 U/L	30 U/L	≤ 3 %
AST STAT	≤ 0.9 U/L	30 U/L	≤ 3 %
AST/P-5-P	≤ 1.2 U/L	40 U/L	≤ 3 %
Bilirubin Total Gen.3	≤ 0.03 mg/dL	0.99 mg/dL	≤ 3 %
Bilirubin, Direct	≤ 0.012 mg/dL	0.30 mg/dL	≤ 3 %
Calcium Gen.2 (serum, urine, STAT, routine)	≤ 0.16 mg/dL	OR	≤ 2 %
Cholesterol (Abell-Kendall)	≤ 3.87 mg/dL	201 mg/dL	≤ 2 %
CK	≤ 3 U/L	140 U/L	≤ 2 %
CO2-L	≤ 0.7 mmol/L	22 mmol/L	≤ 3 %
Creatinine Plus ver 2 (serum)	≤ 0.02 mg/dL	0.90 mg/dL	≤ 2 %
Creatinine Plus ver 2 (urine)	≤ 0.85 mg/dL	28.3 mg/dL	≤ 3 %
Creatinine/Jaffe, rate blanked, compensated (serum)	≤ 0.03 mg/dL	0.90 mg/dL	≤ 3 %
Creatinine/Jaffe, rate blanked (urine)	≤ 0.85 mg/dL	28.3 mg/dL	≤ 3 %
Creatinine/Jaffe STAT (serum)	≤ 0.03 mg/dL	0.90 mg/dL	≤ 3 %
Creatinine/Jaffe STAT (urine)	≤ 0.85 mg/dL	28.3 mg/dL	≤ 3 %
Ethyl Alcohol (serum/plasma/urine)	≤ 1.84 mg/dL	100 mg/dL	≤ 2 %
Fructosamine	≤ 6.0 µmol/L	285 µmol/L	≤ 2 %
GGT ver 2	≤ 0.8 U/L	40 U/L	≤ 2 %

cobas c 701 and c 702 analyzers Within Run Precision Guidelines

Clinical Chemistry

Assay	SD	Sample Concentration	%CV Guideline
Glucose/HK gen. 3 (serum)	≤ 1.44 mg/dL	70.3 mg/dL	≤ 2 %
Glucose/HK gen. 3 STAT (serum)	≤ 1.44 mg/dL	70.3 mg/dL	≤ 2 %
Glucose/HK ver 3 (CSF)	≤ 0.72 mg/dL	39.6 mg/dL	≤ 2 %
Glucose/HK ver 3 STAT (CSF)	≤ 0.72 mg/dL	39.6 mg/dL	≤ 2 %
Glucose/HK ver 3 (urine)	≤ 0.54 mg/dL	19.8 mg/dL	≤ 3 %
Glucose/HK ver 3 STAT (urine)	≤ 0.54 mg/dL	19.8 mg/dL	≤ 3 %
HDL-C gen 4	---	---	≤ 3 %
Iron gen 2	≤ 3.0 µg/dL	150 µg/dL	≤ 2 %
Lactate gen 2 (plasma)	≤ 0.36 mg/dL	19.8 mg/dL	≤ 2 %
Lactate gen 2 STAT (plasma)	≤ 0.36 mg/dL	19.8 mg/dL	≤ 2 %
Lactate gen 2 (CSF)	≤ 0.63 mg/dL	19.8 mg/dL	≤ 3 %
Lactate gen 2 STAT (CSF)	≤ 0.63 mg/dL	19.8 mg/dL	≤ 3 %
LDH ver. 2	≤ 4 U/L	200 U/L	≤ 2 %
LDL-C gen 3	≤ 3.09 mg/dL	155 mg/dL	≤ 2 %
Lipase	≤ 1.2 U/L	60 U/L	≤ 2 %
Lipase STAT	≤ 1.2 U/L	60 U/L	≤ 2 %
Magnesium gen. 2 (serum)	≤ 0.05 mg/dL	1.7 mg/dL	≤ 2 %
Magnesium gen. 2 STAT (serum)	≤ 0.05 mg/dL	1.7 mg/dL	≤ 2 %
Magnesium gen. 2 (urine)	≤ 0.17 mg/dL	4.13 mg/dL	≤ 3 %
Magnesium gen. 2 STAT (urine)	≤ 0.17 mg/dL	4.13 mg/dL	≤ 3 %
Phosphorus (serum)	≤ 0.06 mg/dL	2.70 mg/dL	≤ 2 %
Phosphorus STAT (serum)	≤ 0.06 mg/dL	2.70 mg/dL	≤ 2 %
Phosphorus (urine)	≤ 1.21 mg/dL	40.3 mg/dL	≤ 3 %
Phosphorus STAT (urine)	≤ 1.21 mg/dL	40.3 mg/dL	≤ 3 %
Total Protein	≤ 0.13 g/dL	6.6 g/dL	≤ 2 %
Total STAT Protein	≤ 0.13 g/dL	6.6 g/dL	≤ 2 %
Total Protein U/CSF (CSF)	≤ 1.4 mg/dL	45 mg/dL	≤ 3 %
Total Protein U/CSF (Urine)	≤ 0.36 mg/dL	12.0 mg/dL	≤ 3 %
Triglycerides/GPO	≤ 4.43 mg/dL	204 mg/dL	≤ 2 %
Triglycerides/GB	≤ 4.43 mg/dL	204 mg/dL	≤ 2 %
UIBC (serum/plasma)	≤ 13.4 µg/dL	335 µg/dL	≤ 4 %
Urea/BUN (serum)	≤ 0.48 mg/dL	23.2 mg/dL	≤ 2 %
Urea/BUN STAT (serum)	≤ 0.48 mg/dL	23.2 mg/dL	≤ 2 %
Urea/BUN (urine)	≤ 14 mg/dL	420 mg/dL	≤ 3 %
Urea/BUN STAT (urine)	≤ 14 mg/dL	420 mg/dL	≤ 3 %

cobas c 701 and c 702 analyzers Within Run Precision Guidelines

Clinical Chemistry

Assay	SD	Sample Concentration	%CV Guideline
Uric Acid ver. 2 (serum)	≤ 0.1 mg/dL	7.0 mg/dL	≤ 2 %
Uric Acid ver. 2 (urine)	≤ 2.8 mg/dL	92 mg/dL	≤ 3 %

n = 21 for all tests

cobas c 701 and c 702 analyzers Within Run Precision Guidelines

Ion Selective Electrodes (ISE)

Assay	SD	Sample Concentration	%CV Guideline
Chloride (serum/plasma)	≤ 1.7 mmol/L	100 mmol/L	≤ 1.7%
Chloride (urine)	≤ 2.2 mmol/L	110.0 mmol/L	≤ 2.0%
Potassium (serum/plasma)	≤ 0.05 mmol/L	4.0 mmol/L	≤ 1.2%
Potassium (urine)	≤ 0.5 mmol/L	25 mmol/L	≤ 2.0%
Sodium (serum/plasma)	≤ 1.35 mmol/L	135.0 mmol/L	≤ 1.0%
Sodium (urine)	≤ 1.2 mmol/L	20 - 40 mmol/L	---
	---	40 - 60 mmol/L	≤ 3.0%
		> 60 mmol/L	≤ 2.0%

n = 21 for all tests

cobas c 701 and c 702 analyzers Within Run Precision Guidelines

Specific Proteins

Assay	SD	Sample Concentration	%CV Guideline
α1-Acid Glycoprotein Gen.2	≤ 0.05 g/L	1.2 g/L	≤ 4%
α1-Antitrypsin, ver.2	≤ 4 mg/dL	90 mg/dL	≤ 4%
Antistreptolysin O	≤ 4 IU/mL	100 IU/mL	≤ 4 %
Apolipoprotein A-1 ver.2	≤ 4 mg/dL	100 mg/dL	≤ 4%
Apolipoprotein B ver.2	≤ 4 mg/dL	100 mg/dL	≤ 4%
Albumin TQ (urine)	≤ 0.8 mg/L	20 mg/L	≤ 4%
Albumin TQ (serum)	≤ 1.4 g/L	35 g/L	≤ 4%
Albumin TQ (CSF)	≤ 10 mg/L	240 mg/L	≤ 4%
β2-Microglobulin (serum/plasma)	≤ 0.09 mg/L	2.2 mg/L	≤ 4%
β2-Microglobulin (urine)	≤ 0.05 mg/L	1.0 mg/L	≤ 4%
C-Reactive Protein Gen.3	≤ 0.08 mg/L	0.2 - 1.4 mg/L	---
	---	> 1.4 - 5 mg/L	≤ 5 %
	---	> 5 mg/L	≤ 3 %
Tina-quant C-Reactive Protein IV	---	3 - 5 mg/L	≤ 5 %
	---	> 5 mg/L	≤ 3 %
Cardiac C-Reactive Protein (Latex) High Sensitive	≤ 0.04 mg/L	1.0 mg/L	≤ 4%
Ceruloplasmin	≤ 1 mg/dL	30 mg/dL	≤ 4%
Complement C3c ver.2	≤ 4 mg/dL	90 mg/dL	≤ 4%
Complement C4 ver.2	≤ 0.4 mg/dL	10 mg/dL	≤ 4%
Ferritin Gen.4	< 1.5 µg/L	15 µg/L	---
	---	15 - 40 µg/L	< 10%
	---	> 40 µg/L	< 5%
Haptoglobin ver.2	≤ 1.2 mg/dL	30 mg/dL	≤ 4%
Homocysteine Enzymatic Assay	≤ 0.75 µmol/L	15 µmol/L	≤ 5%
IgA Gen.2 standard	≤ 3.0 mg/dL	70 mg/dL	≤ 4 %
IgA Gen.2 sensitive	≤ 2.0 mg/dL	40 mg/dL	≤ 4 %
IgG Gen.2 standard	≤ 30 mg/dL	700 mg/dL	≤ 4 %
IgG Gen.2 sensitive (CSF)	≤ 0.08 mg/dL	1.5 mg/dL	≤ 5 %
IgM Gen.2 standard	≤ 2.0 mg/dL	40 mg/dL	≤ 4 %
IgM Gen.2 sensitive	≤ 1.0 mg/dL	20 mg/dL	≤ 4 %
Lipoprotein (a) Gen.2	≤ 1.2 mg/dL	30 mg/dL	≤ 4%
Myoglobin Gen.2	≤ 2.4 ng/mL	OR	≤ 4%
Prealbumin	≤ 2 mg/dL	40 mg/dL	≤ 4%
Rheumatoid Factors II	≤ 0.8 IU/mL	14 IU/mL	≤ 6 %
Soluble Transferrin Receptor (sTfR)	≤ 0.1 mg/L	2 mg/L	≤ 4%

cobas c 701 and c 702 analyzers Within Run Precision Guidelines

Specific Proteins

Assay	SD	Sample Concentration	%CV Guideline
Transferrin ver.2 (serum/plasma)	≤ 10 mg/dL	200 mg/dL	≤ 4%
Transferrin ver.2 (urine)	≤ 0.25 mg/L	5 mg/L	≤ 5%

n = 21 for all tests