The background of the page is a warm, orange-toned microscopic image showing cellular structures, including a prominent DNA double helix in the lower-left quadrant and various clusters of cells throughout the field.

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Common genetic alterations	Type I	Type II
PTEN mutation	52~78%	1~11%
PIK3CA mutation	36~52%	24~42%
PIK3R1 mutation	21~43%	0~12%
KRAS mutation	15~43%	2~8%
ARID1A mutation	25~48%	6~11%
CTNNB1 mutation	23~24%	0~3%
TP53 mutation	9~12%	60~91%
HER2 amplification	0	27~44%
Microsatellite instability	28~40%	0~2%

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	POLE	MSI	Copy number low (endometrioid)	Copy number high (serous-like)
Copy number aberrations	Low	Low	Low	High
MSI/MLH1 methylation	Mixed	MSI high	MSI stable	MSI stable
Mutation rate	Very high	High	Low	Low
Genes commonly mutated	POLE (100%)	PTEN (88%)	PTEN (77%)	TP53 (92%)
	PTEN (94%)	RPL22 (37%)	CTNNB1 (52%)	PPP2R1A (22%)
	PIK3CA (71%)	KRAS (35%)	PIK3CA (53%)	PIK3CA (47%)
	PIK3R1 (65%)	PIK3CA (54%)	PIK3R1 (33%)	
	FBXW7 (82%)	PIK3R1 (40%)	ARID1A (42%)	
	ARID1A (76%)	ARID1A (37%)		
	KRAS (53%) ARID5B (47%)			
Histological type	Endometrioid	Endometrioid	Endometrioid	Serous, Endometrioid

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	POLE	MSI	Copy number low (endometrioid)	Copy number high (serous-like)
Histological type	Endometrioid	Endometrioid	Endometrioid	Serous, Endometrioid
Associated histological feature	Ambiguous morphology	Lymphovascular space invasion	Squamous differentiation	Lymphovascular space invasion
	Broad front invasion	MELF-type invasion	ER/PgR expression	Destructive invasion
	TILs, peri-tumoral lymphocytes	TILs, peri-tumor lymphocytes		Giant tumoral cells
	Giant tumoral cells	Low uterine segment involvement		Hobnailing
				Slit-like spaces
Prognosis in early stage (I~II)	Excellent	Intermediate	Excellent/Intermediate/Poor	Poor
Suggested treatment options in recurrent/metastatic disease	Checkpoint inhibitors	Checkpoint inhibitors	Hormonal therapy mTOR inhibitors	Small molecule activators of p53

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- 1) Murali, R., Soslow, R.A., Weigelt, B. : Classification of endometrial carcinoma: more than two types. *Lancet Oncol.* 2014, 15; e268-278
- 2) McAlpine, J., Leon-Castillo, A., Bosse, T. : The rise of novel classification system for endometrial carcinoma; integration of molecular subclasses. *J Pathol.* 2018, 244; 538-549
- 3) Djordjevic, B., Hennessy, B.T., Li, J. et al. : Clinical assessment of PTEN loss in endometrial carcinoma: immunohistochemistry outperforms gene sequencing. *Mod Pathol.* 2012, 25; 699-708
- 4) Kurman, R.J., Carangiu, M.L., Herrington, C.S. et al. (eds.) : *WHO Classification of Tumor of Female Reproductive Organs*, 4th ed., IARC Press, Lyon, 2014
- 5) Mutter, G.L. : Diagnosis of premalignant endometrial disease. *J Clin Pathol.* 2002, 55; 326-331
- 6) Monte, N.M., Webster, K.A., Neuberger, D. et al. : Joint loss of PAX2 and PTEN expression in endometrial precancers and cancer. *Cancer Res*, 2010, 70; 6225-6232
- 7) Huang, E.C., Mutter, G.L., Crum, C.P., Nucci, M.R. : Clinical outcome in diagnostically ambiguous foci of gland crowding in the endometrium. *Mod Pathol* . 2010, 23; 1486-1491
- 8) Mutter, G.L., Ince, T.A., P, Jan. et al. : Molecular identification of latent precancers in histologically normal endometrium. *Cancer Res* 2001, 61; 4311-4314
- 9) Lorito, A.D., Zappacosta, R., Capanna, S. et al. : Expression of PTEN in endometrial liquid-based cytology. *Acta Cytologica* 2014, 58; 495-500

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