


# NTRK Gene Fusion Detection

## Comparison of Available Technologies



	NGS	IHC	FISH	RT-PCR
<b>Ability to detect different fusions</b>	<p><b>Detects fusions in all three NTRK genes<sup>3,6</sup></b> NGS is a sensitive and specific way to detect all the possible variations of NTRK gene fusions and also identify other biomarkers in one comprehensive molecular test.</p>	<p><b>May detect proteins encoded by NTRK genes but requires confirmation by other tests<sup>2,4,6,9</sup></b> To be used as a stratifying tool; other methodologies are required to confirm gene fusions. Pan-TRK antibodies may detect presence of protein encoded by any of the three NTRK genes but can not differentiate between wild type and fusion proteins<sup>2,4,9</sup></p>	<p><b>Can detect NTRK fusions, depending on the probes used<sup>11,12</sup></b> Different probes are required for each NTRK gene. For ROS1 fusion detection dual colour 'break-apart' probes are used which label the 3' (centromeric) part of the fusion breakpoint with one colour fluorochrome and the 5' (telomeric) part with another. Rearrangements are determined by looking at the pattern of fluorochrome expression.<sup>6,13</sup></p>	<p><b>Requires multiple primer sets to detect gene fusions<sup>11,14</sup></b> As the location of gene rearrangements is not known at the time of testing, multiple primer sets are required, one for each fusion variant.<sup>11,13,14</sup></p>
<b>Reliability</b>	<p><b>Reliability varies depending on the assay used<sup>10</sup></b> Different NGS panels target different regions of the sequence and depth of coverage varies between assays. Requirement for tumor content can also vary between assays.<sup>10</sup></p>	<p><b>Detects NTRK proteins with high specificity and sensitivity<sup>2,4</sup></b> Pan-TRK antibodies have been reported to have 95–100% sensitivity and up to 100% specificity for TRK proteins depending on the threshold used to define positivity.<sup>2,6,7</sup></p>	<p><b>False positives or negatives might occur<sup>6,11</sup></b> Reliability depends on the probes used. There is a risk of false-positive results due to complex chromosomal translocations and detection of non-functional fusion proteins.<sup>6,11</sup> Variant or complex rearrangements may be missed.<sup>1,2,6,11</sup> False negatives results may be above 30% in some cases.<sup>11</sup></p>	<p><b>Reliable for known gene fusions, requires FISH confirmation when testing for NTRK fusions<sup>7</sup></b> Each variant requires a specific primer set. Therefore, PCR may miss unknown or untested variants.<sup>5,6,13,14</sup> ROS1 expression is detectable in normal lung tissue, PCR testing is not specific enough to detect ROS1-positive lung cancer without FISH confirmation.<sup>8,14</sup></p>
<b>Advantage</b>	<p><b>Can detect novel fusion partners and evaluates multiple actionable targets<sup>11</sup></b> Depending on the assay used, NGS can evaluate multiple and novel fusion partners while preserving limited tissue. RNA sequencing is focused on coding sequences rather than introns, and is suited to gene fusion detection.<sup>4,11</sup></p>	<p><b>Low cost and readily available<sup>2,6,11</sup></b> Pan-TRK antibodies detect TRK A, B and C proteins with a turnaround time of 1–2 days.<sup>11</sup></p>	<p><b>Readily available</b> Allows visualisation of the target within the cell and enables several targets to be detected in one sample using multiple fluorophores.<sup>7,11</sup> The use of break-apart probes allows fusions with unknown partners to be detected.<sup>11</sup></p>	<p><b>Low cost per assay with high sensitivity and specificity<sup>11</sup></b></p>
<b>Disadvantage</b>	<p><b>May not identify all NTRK gene fusions<sup>11</sup></b> DNA-NGS is limited by intron size. RNA-NGS is suited to gene fusion detection but may be limited by RNA quality.<sup>4,11</sup></p>	<p><b>Cannot reliably differentiate between normal protein expression and proteins resulting from gene rearrangements<sup>6,11</sup></b> IHC can detect both wild-type and fusion TRK proteins leading to possible false positives. There is also a risk of possible false negatives for fusions involving TRKC. Therefore, orthogonal techniques are required to confirm the presence of gene rearrangements.<sup>6,11</sup></p>	<p><b>Conventional FISH may require multiple tests to detect NTRK fusions</b> The target sequence must be known for conventional FISH and three separate tests are required for NTRK1, NTRK2 and NTRK3.<sup>11</sup> Depending on the type of translocation there is also a risk of false positive or false negative results.<sup>6,11</sup></p>	<p><b>Can only detect known target sequences<sup>11</sup></b> Target sequences must be known and cannot readily detect novel fusion partners. A comprehensive multiplex reverse transcriptase polymerase chain reaction (RT-PCR) assay might be challenging because of the potentially large number of possible 5' fusion partners.<sup>7,11</sup></p>

## References

1. Vaishnavi A, Le AT, Doebele RC. TRKING down an old oncogene in a new era of targeted therapy. *Cancer Discov.* 2015;5(1):25-34.
2. Murphy DA, Ely HA, Shoemaker R, et al. Detecting gene rearrangements in patient populations through a 2-step diagnostic test comprised of rapid IHC enrichment followed by sensitive next-generation sequencing. *Appl Immunohistochem Mol Morphol.* 2017;25(7):513-523.
3. Su D, Zhang D, Chen K, et al. High performance of targeted next generation sequencing on variance detection in clinical tumor specimens in comparison with current conventional methods. *J Exp Clin Cancer Res.* 2017;36(1):121.
4. Du X, Shao Y, Qin HF, Tai YH, Gao HJ. ALK-rearrangement in non-smallcell lung cancer (NSCLC). *Thorac Cancer.* 2018;9(4):423-430.
5. Ali G, Bruno R, Savino M, et al. Analysis of fusion genes by NanoString system: a role in lung cytology? *Arch Pathol Lab Med.* 2018;142(4):480-489.
6. Precision Oncology News. Japanese Regulator Approves Foundation Medicine CDx Assay for Roche's Rozlytrek. <https://www.precisiononcologynews.com/regulatory-news/japaneseregulator-approves-foundation-medicine-cdx-assay-roches-rozlytrek#.Xj2faBd7ITY>. June 27, 2019. Accessed February 3, 2020.
7. Penault-Llorca F, Rudzinski ER, Sepulveda AR. Testing algorithm for identification of patients with TRK fusion cancer. *J Clin Pathol.* 2019;72(7):460-467.
8. Referenced with permission from the NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) for Non-Small Cell Lung Cancer. V.2.2020. © 2019 National Comprehensive Cancer Network, Inc. All rights reserved. Published December 23, 2019. Accessed February 7, 2020. To view the most recent and complete version of the guideline, go online to [www.NCCN.org](http://www.NCCN.org).
9. Roche Tissue Diagnostics. VENTANA pan-TRK (EPR17341) Assay: Package Insert 1017533EN Rev A. [https://productlibrary.ventana.com/ventana\\_portal/executeSearch.do](https://productlibrary.ventana.com/ventana_portal/executeSearch.do). Accessed February 7, 2020.
10. Wang Z, Gerstein M, Snyder M. RNA-Seq: a revolutionary tool for transcriptomics. *Nat Rev Genet.* 2009;10(1):57-63. doi:10.1038/nrg2484.
11. Aguado C, Giménez-Capitán A, Karachaliou N, et al. Fusion gene and splice variant analyses in liquid biopsies of lung cancer patients. *Transl Lung Cancer Res.* 2016;5(5):525-531.
12. Foundation Medicine. FoundationOne® Liquid: Technical Specifications. <https://www.foundationmedicine.com/genomic-testing/foundation-one-liquid>. Accessed February 7, 2020.
13. Foundation Medicine. FoundationOne® Heme: Technical Specifications. <https://www.foundationmedicine.com/genomic-testing/foundation-one-heme>. Accessed February 7, 2020.
14. Burd EM. Validation of laboratory-developed molecular assays for infectious diseases. *Clin Microbiol Rev.* 2010;23(3):550-576.

To learn more about pan-TRK IHC testing, visit [go.roche.com/NTRKIHC](https://go.roche.com/NTRKIHC) or contact your Roche representative.

All product names and trademarks are the property of their respective owners.

© 2020 Roche. MC-US-07486-1020

Roche Diagnostics Corporation  
9115 Hague Road  
Indianapolis, Indiana 46256

[diagnostics.roche.com](https://diagnostics.roche.com)