Urine Specimen Collection Guide

Read the collection instructions. If you have any questions about how to collect the sample, please ask your healthcare provider.

Preparing for sample collection

- Wash hands prior to collection.

Please follow these steps to assure reliable results and reduce the need for retesting.

Note

- Patient should not urinate 1 hour prior to sampling.
- It may be helpful to mark the volume required on the cup.
- For best results, female patients should not cleanse the labial area prior to collection.

Collecting male/female urine

- 

Prepare
Position the urine cup to collect the beginning of the urine stream.
- 


Collect urine
Begin urinating. Collect the first 10–50 ml of urine.
- 

Secure container
Close the lid of the urine cup securely, label and date.


How to transfer urine specimen to cobas® PCR tube

- 

Collect
Prior to sampling, the patient should not have urinated for at least one hour. Given that collection of larger volumes of urine may reduce test sensitivity, please direct patient to provide first-catch urine (approximately 10 to 50 mL of the initial urine stream) into a urine collection cup (not provided).
- 

Pipette
Immediately transfer urine into the **cobas®** PCR Media tube using the provided disposable pipette.
- 

Transfer
The correct volume of urine has been added when the fluid level is between the two black lines on the tube label.
- 

Cap
Tightly re-cap the **cobas®** PCR Media tube.
- 

Mix
Invert the tube 5 times to mix. The specimen is now ready for transport.

How to transport samples for testing

Urine collection Kit

Handling precaution

The collection tube media can cause irritation if contacted with skin or other body parts. Handle the collection tube carefully.

- Use care to avoid splashing contents of the tube. If the contents of the tube are spilled on your skin, wash the affected area with soap and water. If the contents of the tube splash into your eyes, flush them with water immediately. Always notify your healthcare provider.

Note

In case the contents of the tube are accidentally spilled, do not attempt to clean up. Immediately notify your healthcare provider for appropriate action. If the urine specimen cannot be transferred immediately to the **cobas** PCR media tube, it can be stored at 2°C to 30°C for up to 24 hours.

cobas® PCR Media tube: Shipping and storage

cobas x800 users: Urine specimens collected in **cobas®** PCR Media are stable at 2–30° C for up to 1 year.

Liat users: Urine samples may be stored at 2–30°C for up to 3 hours after collection, or at 2–8°C for up to 48 hours after collection.

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