


# **VENTANA PD-L1 (SP142) Assay**

## *Guiding immunotherapy in NSCLC*

**PD-L1**




Hiker's path: VENTANA PD-L1 (SP142) Assay  
on non-small cell lung cancer tissue  
Location: Point Conception, CA



# VENTANA PD-L1 (SP142) Assay

*Assess NSCLC patient benefit from TECENTRIQ*



Positive NSCLC tissue stained with PD-L1 (SP142) assay, 10x

## Intended use statement

VENTANA PD-L1 (SP142) Assay is a qualitative immunohistochemical assay using rabbit monoclonal anti-PD-L1 clone SP142 intended for use in the assessment of the PD-L1 protein in formalin-fixed, paraffin-embedded (FFPE) urothelial carcinoma and non-small cell lung cancer (NSCLC) tissue stained with OptiView DAB IHC Detection Kit and OptiView Amplification Kit on a VENTANA BenchMark ULTRA instrument. Determination of PD-L1 status is indication-specific, and evaluation is based on either the proportion of tumor area occupied by PD-L1 expressing tumor-infiltrating immune cells (% IC) of any intensity or the percentage of PD-L1 expressing tumor cells (% TC) of any intensity.

PD-L1 expression in  $\geq 5\%$  IC determined by VENTANA PD-L1 (SP142) Assay in urothelial carcinoma tissue is associated with increased objective response rate (ORR) in a non-randomized study of TECENTRIQ® (atezolizumab).

PD-L1 expression in  $\geq 50\%$  TC or  $\geq 10\%$  IC determined by VENTANA PD-L1 (SP142) Assay in NSCLC tissue may be associated with enhanced overall survival from TECENTRIQ (atezolizumab).

## PD-L1 diagnostic confidence


Using the right test to determine PD-L1 status for immunotherapy options is important, and the VENTANA PD-L1 (SP142) Assay is the only FDA approved test for TECENTRIQ. This innovative assay is the first to evaluate patient PD-L1 expression using both tumor cell (TC) and tumor-infiltrating immune cell (IC) staining. Determining a patient's PD-L1 expression level can give insight to the overall survival that may be achieved from TECENTRIQ.\*

The VENTANA PD-L1 (SP142) Assay

- FDA approved to assess NSCLC patient treatment benefit from TECENTRIQ
- Informative for the clinician of a patient's potential overall survival
- Novel scoring algorithm using PD-L1 staining in both TC and IC
- Designed to enhance visual contrast of immune cell staining within the tumor microenvironment

The PD-L1 (SP142) Assay gives you the confidence to guide immunotherapy decisions in NSCLC.

*\*All randomized patients in a NSCLC phase III study observed benefit from TECENTRIQ regardless of PD-L1 status.*



**Rabbit**  
**Monoclonal**

**IHC antibody developed by**  
**Spring Bioscience Corporation**



**Fully**  
**Automated**

**With specific and**  
**robust signal**



**Highly**  
**Reproducible**

**Accurate scoring and**  
**educational resources**

## About PD-L1

PD-L1 is a transmembrane protein that down-regulates immune responses through binding to its two inhibitory receptors, programmed death-1 (PD-1) and B7.1. PD-1 is an inhibitory receptor expressed on T cells following T-cell activation, which is sustained in states of chronic stimulation such as in chronic infection or cancer.<sup>1</sup> Ligation of PD-L1 with PD-1 inhibits T-cell proliferation, cytokine production and cytolytic activity, leading to the functional inactivation or exhaustion of T cells. B7.1 is a molecule expressed on antigen presenting cells and activated T cells. PD-L1 binding to B7.1 on T cells and antigen presenting cells can mediate down-regulation of immune responses, including inhibition of T-cell activation and cytokine production.<sup>2</sup> PD-L1 expression has been observed in immune cells and tumor cells.<sup>3,4</sup> Aberrant expression of PD-L1 on tumor cells has been reported to impede anti-tumor immunity, resulting in immune evasion.<sup>1</sup> Therefore, interruption of the PD-L1/PD-1 pathway represents an attractive strategy to reinvigorate tumor-specific T-cell immunity suppressed by the expression of PD-L1 in the tumor microenvironment.


## About non-small cell lung cancer

Lung cancer has been the most common cancer in the world for several decades and remains the leading cause of cancer deaths worldwide. It is estimated to account for 12.9% of all new cancer cases and is responsible for nearly 1.59 million deaths annually worldwide, or approximately one in five cancer-related deaths.<sup>5</sup> In the United States alone, approximately 224,390 new cases of lung cancer and 158,080 lung cancer-related deaths are predicted for 2016.<sup>6</sup> Although improvements have been made in diagnosis and therapy options, prognosis remains poor with low long-term survival rates for all stages. Over the past three decades, lung cancer has shown among the least improvement in survival rates when compared with other cancers.<sup>7</sup>

Non-small cell lung cancer (NSCLC), one of the two major types of lung cancer, accounts for approximately 85% of all lung cancer cases.<sup>8</sup> In more than half of patients newly diagnosed with NSCLC, the disease has already metastasized, greatly decreasing the likelihood of survival. The five-year relative survival rate for NSCLC diagnosed as distant disease is 4.7%.<sup>8</sup>


The majority of patients with NSCLC present with inoperable, locally advanced disease (Stage IIIB) or metastatic disease (Stage IV), neither of which currently has any curative treatment options. On average, these patients die within a year of diagnosis. Improvement in the clinical outcome of lung cancer is likely to be achieved through better understanding of the molecular events that underlie its pathogenesis, identifying new biomarker targets and developing new treatment options.

## The PD-L1 immunologic checkpoint



## The tumor microenvironment

*A key player in tumor progression*



*The PD-L1 (SP142) Assay stain highlights a heterogeneous population of immune cells. The majority of these cells are morphologically consistent with lymphocytes, macrophages, dendritic cells and granulocytes.*

### *What is the tumor microenvironment (TME)?*

The tumor microenvironment (TME) consists of different cells, including immune cells and non-cellular components in and around the tumor. The TME has been recognized to play a significant role in tumor progression.<sup>9,10</sup>

### *Why is the TME important?*

The TME shapes tumor evolution (whether the tumor regresses, develops resistance, evades the immune system and/or metastasizes) and consequently impacts patient outcomes.<sup>11</sup> An association has been observed between the levels of tumor infiltrating immune cells, key components of the TME and outcome in lung cancer patients. A meta-analysis showed that high levels of CD8+ lymphocyte tumor stroma infiltration was associated with better overall survival (Hazard ratio = 0.76 95% confidence interval 0.62-0.93), high levels of CD4+ lymphocyte tumor stroma infiltration was also associated with better overall survival (Hazard ratio = 0.65 95% confidence interval 0.46-0.91, P 0.013). In contrast, high levels of FOXP3P+ T lymphocyte stroma infiltration was associated with poor progression-free survival (Hazard ratio = 2.67 95% confidence interval 1.74-4.08, P value < 0.001).<sup>12</sup>

### *What is the role of PD-L1 in the TME?*

Aberrant expression of PD-L1 on tumor cells has been reported to impede anti-tumor immunity, resulting in immune evasion.<sup>13</sup> Therefore, interruption of the PD-L1/PD-1 pathway represents an attractive strategy to reinvigorate tumor-specific T-cell immunity suppressed by the expression of PD-L1 in the tumor and in the tumor microenvironment. This approach has proven effective: PD-L1 expression on immune cells in the tumor microenvironment has been shown to help assess NSCLC patient benefit from TECENTRIQ, an anti-PD-L1 drug.<sup>14</sup>





# Scoring of VENTANA PD-L1 (SP142) Assay

## General guidelines



Immunohistochemical (IHC) staining with the VENTANA PD-L1 (SP142) Assay demonstrates staining in tumor cells (TC, Figure 1) as well as tumor-infiltrating immune cells (IC, Figure 2).

A patient’s tumor tissue stained with the VENTANA PD-L1 (SP142) Assay is evaluated for TC and/or IC staining. TC are scored as the proportion of viable tumor cells showing PD-L1 membrane staining of any intensity. IC are scored as the proportion of tumor area that is occupied by PD-L1 staining IC of any intensity.




## Differentiation of TC and IC Staining

TC staining can be observed in association with IC staining. Review of the corresponding H&E slide will help in identifying IC among TC. This along with a high magnification review of the PD-L1 stained slide may aid in differentiating between TC and IC staining. The following images demonstrate different commonly observed patterns encountered in clinical practice, when TC and IC staining is observed together.



VENTANA PD-L1 (SP142) Assay scoring approach – NSCLC

NSCLC tissue stained with the VENTANA PD-L1 (SP142) Assay will be scored using a stepwise approach according to the criteria outlined in Table 1. TC are scored as the proportion of viable tumor cells showing PD-L1 membrane staining of any intensity. IC are scored as the proportion of tumor area that is occupied by PD-L1 staining IC of any intensity. VENTANA PD-L1 (SP142) Assay-stained slides will first be evaluated for TC staining (Step 1 in Table 1). If the specimen contains any discernible PD-L1 membrane staining of any intensity in  $\geq 50\%$  TC, the case will be assigned a PD-L1 expression level of  $\geq 50\%$  TC. If the specimen contains  $< 50\%$  TC staining, the slide will then be evaluated for IC staining (Step 2 in Table 1). If the specimen contains PD-L1 staining of any intensity in IC occupying  $\geq 10\%$  of tumor area, the case will be assigned a PD-L1 expression level of  $> 10\%$  IC. If the specimen contains PD-L1 staining of any intensity in IC covering  $< 10\%$  of tumor area, the case will be assigned a PD-L1 expression level of  $< 50\%$  TC and  $< 10\%$  IC.




Stepwise scoring process – VENTANA PD-L1 (SP142) Assay in NSCLC tissues (images 10x magnification).

Table 1: VENTANA PD-L1 (SP142) Assay Stepwise Scoring Algorithm for NSCLC	
Step 1: Tumor Cell (TC) Staining Assessment	PD-L1 Expression
Presence of discernible PD-L1 membrane staining of any intensity in $\geq 50\%$ of tumor cells	$\geq 50\%$ TC
Absence of any discernible PD-L1 staining -or- Presence of discernible PD-L1 membrane staining of any intensity in $< 50\%$ of tumor cells	Proceed to Step 2
Step 2: Tumor Infiltrating Immune Cell (IC) Staining Assessment	PD-L1 Expression
Presence of discernible PD-L1 staining of any intensity in tumor-infiltrating immune cells covering $\geq 10\%$ of tumor area occupied by tumor cells, associated intratumoral and contiguous peritumoral stroma	$\geq 10\%$ TC
Absence of any discernible PD-L1 staining -or- Presence of discernible PD-L1 staining of any intensity in tumor-infiltrating immune cells covering $< 10\%$ of tumor area occupied by tumor cells, associated intratumoral and contiguous peritumoral stroma	$< 50\%$ TC and $< 10\%$ IC

## VENTANA PD-L1 (SP142) Assay results from clinical trials

Results from the OAK trial<sup>\* 14</sup>



- All randomized patients in a NSCLC phase III study observed benefit from TECENTRIQ regardless of PD-L1 status.

Efficacy Results in the Primary Analysis Population from the Phase III Study		
Overall Survival	TECENTRIQ (n=425)	Docetaxel (n=425)
Deaths (%)	271 (64%)	298 (70%)
Median, months (95% CI)	13.8 (11.8, 15.7)	9.6 (8.6, 11.2)
Hazard ratio* (95% CI)	0.74 (0.63, 0.87)	
p-value**	0.0004	

\*Stratified by PD-L1 expression in tumor-infiltrating immune cells, the number of prior chemotherapy regimens, and histology

\*\*Based on the stratified log-rank test

CI = confidence interval

Tumor specimens were evaluated prospectively using VENTANA PD-L1 (SP142) Assay at a central laboratory and the results were used to define the PD-L1 expression subgroups for pre-specified analyses. Of the 850 patients, 16% were classified as having high PD-L1 expression, defined as having PD-L1 expression on  $\geq 50\%$  of TC or  $\geq 10\%$  of IC. In an exploratory efficacy subgroup analysis of OS based on PD-L1 expression, the hazard ratio was 0.41 (95% CI: 0.27, 0.64) in the high PD-L1 expression subgroup and 0.82 (95% CI: 0.68, 0.98) in patients who did not have high PD-L1 expression.

\*A phase III, open-label, multicenter, randomized study to investigate the efficacy and safety of atezolizumab compared with docetaxel in patients with non-small cell lung cancer after failure with platinum containing chemotherapy

Clinicaltrials.gov number NCT0202008227

	<b>VENTANA PD-L1 (SP142) Assay</b>	<b>OptiView DAB IHC Detection Kit</b>	<b>OptiView Amplification Kit</b>	<b>Negative Reagent Control</b>
<b>Catalog number</b>	<b>740-4859</b> 07709374001	<b>760-700</b> 06396500001	<b>860-099</b> 06718663001	<b>790-4795</b> 06683380001
<b>Quantity</b>	50 tests	250 tests	250 tests	50 tests
<b>Positive control</b>	Tonsil			
<b>Species</b>	Rabbit			
<b>Localization</b>	Membranous and/or Cytoplasmic			

## References

- Blank, C and Mackensen, A, Contribution of the PD-L1/PD-1 pathway to T-cell exhaustion: an update on implications for chronic infections and tumor evasion. *Cancer Immunol Immunother*, 2007. 56(5): p. 739-745.
- Butte MJ, Keir ME, Phamduy TB, et al. Programmed death-1 ligand 1 interacts specifically with the B7-1 costimulatory molecule to inhibit T cell responses. *Immunity*. 2007;27(1):111-122.
- Dong H, Zhu G, Tamada K, Chen L. B7-H1, a third member of the B7 family, co-stimulates T-cell proliferation and interleukin-10 secretion. *Nat Med*. 1999;5(12):1365-1369.
- Herbst RS, Soria JC, Kowanetz M, et al. Predictive correlates of response to the anti-PD-L1 antibody MPDL3280A in cancer patients. *Nature*. 2014;515(7528):563-567.
- Ferlay J, Soerjomataram I, Ervik M, et al. GLOBOCAN 2012 v1.0. Cancer Incidence and Mortality Worldwide: IARC CancerBase No. 11. <http://globocan.iarc.fr>. Published 2013-12-12. Updated 2014-01-09. Accessed 2016-02-08.
- American Cancer Society. *Cancer Facts & Figures* 2016. Atlanta, GA: American Cancer Society; 2016.
- Siegel R, Naishadham D, Jemal A. *Cancer Statistics*, 2013. *CA Cancer J Clin*. 2013;63(1):11-30.
- Howlander N, Noone AM, Krapcho M, et al. (eds). *SEER Cancer Statistics Review (CSR), 1975-2012*. National Cancer Institute. [http://seer.cancer.gov/csr/1975\\_2012/](http://seer.cancer.gov/csr/1975_2012/). Published 2015-04-23. Updated 2015-11-18. Accessed 2016-02-08.
- NCI Dictionary of Cancer Terms. Tumor microenvironment. <http://www.cancer.gov/publications/dictionaries/cancer-terms?cdrid=561725>
- AACR. The Function of Tumor Microenvironment in Cancer Progression. <http://www.aacr.org/Meetings/Pages/MeetingDetail.aspx?EventItemID=73#>. V6pCFPrKaE
- Chen F, Zhuang X, Lin L, Yu P, Wang Y, Shi Y, Hu G, Sun Y. *BMC Med*. New horizons in tumor microenvironment biology: challenges and opportunities. 2015 Mar 5;13:45. doi: 10.1186/s12916-015-0278-7.
- Geng Y, Shao Y, He W, Hu W, Xu Y, Chen J, Wu C, Jiang J. Prognostic Role of Tumor-Infiltrating Lymphocytes in Lung Cancer: a Meta-Analysis. *Cell Physiol Biochem*. 2015;37(4):1560-71. doi: 10.1159/000438523. Epub 2015 Oct 30.
- Blank, C and Mackensen, A, Contribution of the PD-L1/PD-1 pathway to T-cell exhaustion: an update on implications for chronic infections and tumor evasion. *Cancer Immunol Immunother*, 2007. 56(5): p. 739-745).
- Ventana Medical Systems Inc. VENTANA PD-L1 (SP142) Assay Package Insert.

[usdiagnostics.roche.com](http://usdiagnostics.roche.com)  
[www.ventana.com](http://www.ventana.com)

© 2016 Roche PP-US-09136-1016

Roche Diagnostics Corporation  
9115 Hague Road  
Indianapolis, Indiana 46256

VENTANA, BENCHMARK, OPTIVIEW are trademarks of Roche.  
All other product names and trademarks are the property of their respective owners.