

Establishing Previous Test Result as Current Results


After having repeated a test, with or without dilution, you realize that the last result you obtained is not valid. Therefore, you want to use a previous result as the current one.

To Establish Previous Test Results as Current Results

1. Once logged into the **cobas[®] infinity** software, select the **Main** tab
2. Click on **Menu**
3. Navigate to **Work areas**
4. Select the Work area you wish to view

Please note: The Work area names are customized per laboratory. What you see in our example won't necessarily reflect the names of your Work Areas.

5. Double click the desired order from the Monitoring screen to bring up the Validation screen. When accessing a Validation screen from a work area, consider the following:
 - All the filters selected in the work area are also applied in the Validation screen. Therefore, if work area tests are filtered by location, the same filter is applied in the Validation screen
 - By double-clicking an order in the work area, you apply the criterion selected in the Validation criterion drop-down list. Therefore, the orders matching the criterion are validated



Priority	Entry Date	SampleID	Entry Time	First name	Last name	Patient ID1	WA Status
<input checked="" type="checkbox"/>	01/28/2019	Z40000001	10:38:51	Theo	Edwards	123456	■ ■ ■ ■ ■ ■ ■
<input type="checkbox"/>	06/21/2019	Z00000011	19:41:42	Brennan	Mathews	P123477	■ ■ ■ ■ ■ ■ ■
<input type="checkbox"/>	07/02/2019	Z00000012	09:30:41	Brennan	Mathews	P123477	■ ■ ■ ■ ■ ■ ■
<input type="checkbox"/>	07/02/2019	Z00000014	09:39:01	Estrella	Douglas	P123480	■ ■ ■ ■ ■ ■ ■
<input type="checkbox"/>	07/11/2019	Z00000015	09:31:54	Vern	Smith	000008	■ ■ ■ ■ ■ ■ ■
<input type="checkbox"/>	07/11/2019	Z00000016	10:54:45	Adam	Smith	000009	■ ■ ■ ■ ■ ■ ■

Establishing Previous Test Result as Current Results

6. Select the check box of a test
7. Click the **Repetitions** button

<input type="checkbox"/>				Test name	Test result	Primary unit	T...	Ref. ranges	Critical range	Result alarms	Status	Test alarms	Last repeat...	Instrument
<input type="checkbox"/>				Bicarbonate	23	mmol/L		22 - 29				0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>				BUN	24	mg/dl		8 - 25				0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>				Calcium	7.5	mg/dL		8.0 - 11.0		▼		0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>				Hemolytic	3							0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>				Icteric	27							0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>				Lactate	4.20	mmol/l		0.50 - 2.20	> 4.00	▲▲		0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>				Lipemic	11							0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>				Procalcitonin	2.50	ng/mL		0.00 - 0.08	> 1.00	▲▲		0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>				GLUCOSE	138	mg/dL		65.00 - 99.00	< 30.00, > 500.0	▲		0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>				CREATININE	1	mg/dl		0.70 - 1.30	> 10.00			0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>				SODIUM	142			135.00 - 145.00	< 125.00, > 155.0			0		L3_C8K-1 MU1-e602-1-2
<input checked="" type="checkbox"/>				POTASSIUM	5.80			3.50 - 5.40	< 2.50, > 6.50	▲		0	5.70	L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>				CHLORIDE	101			97.00 - 107.00	< 90.00, > 120.0			0		L3_C8K-1 MU1-e602-1-2

13 records

6

7

Repetitions

8. In the Repetitions screen, select the result you want to establish as the current one
9. Click the **Current** button
10. Return to your Validation screen
11. Click on button
12. The result you selected as current is now displayed as the test result

Overview 8 1 General Lab ROCHE ROCHE 10 Log off

Main Administration Monitoring Search

Menu Work areas: Serum Work Area Validation Repeated results

Order ID	Test ID	Test abbrevia...	Test name	Result date	Result time	Result value	Notifications	Flags	C	Current
Z4000001	107	K	POTASSIUM	01/28/2019	10:58:05	5.80	▲	0		<input checked="" type="checkbox"/>
Z4000001	107	K	POTASSIUM	01/28/2019	10:55:13	5.70	▲	0		<input type="checkbox"/>

page 1 of 1 Records: 2

8

9

Current