


Repeating Tests

1. Once logged into the **cobas[®] infinity** software, select the **Main** tab
2. Click on **Menu**
3. Navigate to **Work areas**
4. Select the Work area you wish to view

Please note: The Work Area names are customized per laboratory. What you see in our example won't necessarily reflect the names of your Work areas. The above navigation will help you see the names of all Work areas you have created.


5. Double click the desired order from the Monitoring screen to bring up the Validation screen. When accessing a Validation screen from a work area, consider the following:
 - All the filters selected in the work area are also applied in the Validation screen. Therefore, if work area tests are filtered by location, the same filter is applied in the Validation screen
 - By double-clicking an order in the work area, you apply the criterion selected in the Validation criterion drop-down list. Therefore, the orders matching the criterion are validated



Priority	Entry Date	SampleID	Entry Time	First name	Last name	Patient ID1	WA Status
	01/28/2019	Z40000001	10:38:51	Theo	Edwards	123456	■■■■■■
	06/21/2019	Z00000011	19:41:42	Brennan	Matthews	P123477	■□□□□□
	07/02/2019	Z00000012	09:30:41	Brennan	Matthews	P123477	■□□□□□
	07/02/2019	Z00000014	09:39:01	Estrella	Douglas	P123480	■□□□□□
	07/11/2019	Z00000015	09:31:54	Vern	Smith	000008	■□□□□□
	07/11/2019	Z00000016	10:54:45	Adam	Smith	000009	■□□□□□

Repeating Tests

6. Select the check box of the test that you want to repeat
7. Click on the  button
8. The **Last repeated res.** column displays the previous test result
9. The **Test result** column displays the new result, when available



Order ID: Z4000001 | 123456 - Edwards, Theo - 57 - Male - 10/05/1961

<input type="checkbox"/>		Test name	Test result	Primary unit	T	Ref. ranges	Critical range	Result alarms	Status	Test alarms	Last repeat...	Instrument
<input type="checkbox"/>		Bicarbonate	23	mmol/L		22 - 29				0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>		BUN	24	mg/dl		8 - 25				0		L3_C8K-1 MU1-e602-1-2
<input checked="" type="checkbox"/>		Calcium	7.5	mg/dL		8.0 - 11.0				0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>		Hemolytic	3							0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>		Icteric	27							0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>		Lactate	4.20	mmol/l		0.50 - 2.20	> 4.00	▲▲		0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>		Lipemic	11							0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>		Procalcitonin	2.50	ng/mL		0.00 - 0.08	> 1.00	▲▲		0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>		GLUCOSE	138	mg/dL		65.00 - 99.00	< 30.00, > 500.00	▲		0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>		CREATININE	1	mg/dl		0.70 - 1.30	> 10.00			0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>		SODIUM	142			135.00 - 145.00	< 125.00, > 155.00			0		L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>		POTASSIUM	5.80			3.50 - 5.40	< 2.50, > 6.50	▲		0	5.70	L3_C8K-1 MU1-e602-1-2
<input type="checkbox"/>		CHLORIDE	101			97.00 - 107.00	< 90.00, > 120.00			0		L3_C8K-1 MU1-e602-1-2

13 records

Repetitions