

## New firmware version 01.09.00 for your CoaguChek INRange / CoaguChek Vantus meter

---

**Date:**

Jan 2022

Rev. 1.0

---

We herewith introduce the new firmware version 01.09.00 for your CoaguChek INRange / CoaguChek Vantus meter.

This software version cannot be downgraded.

The differences between the former firmware version 01.06.02 and the new firmware version 01.09.00 are as follows:

- The display of the battery power level, which is shown with the segmented battery icon on the main menu, was improved. Due to a more precise voltage measurement of the HW variants, it gives now a more accurate status of the battery power level.
- When using the embedded reports, the trend graph range is now also available for a time range of 730 days.

Should you have any question about the CoaguChek INRange / CoaguChek Vantus firmware version 01.09.00, please contact your Roche Service organization.

Please note:

The CoaguChek Vantus is available for the US market only.

Sincerely