

The LOINC codes provided in this document are Roche's best interpretation of the LOINC codes that could be associated with each of our reagents. To our knowledge, the codes provided are accurate and reflect a correct LOINC association. However, we are providing LOINC codes for informational use only. **Please note that Roche does not create or manage LOINC coding assignments.** LOINC codes are created and managed by the Regenstrief Institute, Inc. ("Regenstrief"). Roche cooperates with Regenstrief when asked about unique markers and conducts periodic method reviews to determine whether new LOINC codes should be created for unique markers or markers not properly characterized by existing codes. Moreover, the development, assignment, and management of LOINC codes remains solely within the purview of Regenstrief. With respect to certain tests, there may be multiple LOINC codes that could apply. You must determine the most appropriate code for your organization's test catalog number based on specimen type, challenge, and timing, etc.

**Roche makes no representation and/or warranty with respect to the provision of LOINC codes, including but not limited to the suitability of use to which you intend to put them. Roche further disclaims any and all liability associated with your decision to rely upon the LOINC codes provided herein.** You, as a testing laboratory, remain solely responsible for ensuring the accurate selection of LOINC codes and the accurate mapping of LOINC codes within your laboratory information system.

Please direct any questions about the information in this document to the Roche Support Network Customer Support Center at 1-800-800-5973.

**NOTICE:** The information in this document includes all or a portion of the LOINC table, LOINC panels, and forms file and/or LOINC hierarchy file, or is derived from one or more of the foregoing, subject to a license from Regenstrief Institute, Inc. Your use of the LOINC table, LOINC codes, LOINC panels and forms file, and LOINC hierarchies file is also subject to this license, a copy of which is available at <http://loinc.org/terms-of-use>. The current complete LOINC table, LOINC Users' Guide, LOINC panels and forms file, and LOINC hierarchies file are available for download at <http://loinc.org>. The LOINC table and LOINC codes are copyright © 1995-2017, Regenstrief Institute, Inc. All rights reserved. THE LOINC TABLE (IN ALL FORMATS), LOINC PANELS AND FORMS FILE, AND LOINC HIERARCHIES ARE PROVIDED "AS IS." ANY EXPRESS OR IMPLIED WARRANTIES ARE DISCLAIMED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.



©2025, Roche Diagnostics. All rights reserved.

Distribution in USA by: Roche Diagnostics, Indianapolis, IN.

This document is available on the Roche Diagnostics website at [navifyportal.roche.com](http://navifyportal.roche.com).

LOINC is a registered United States trademark of Regenstrief Institute, Inc.

COBAS, AMPLIPREP, TAQMAN, LIGHTCYCLER, LIAT, and MRSA ADVANCED are trademarks of Roche.

The technology used for LightCycler is licensed from Idaho Technology, Inc.

All other product names and trademarks are the property of their respective owners.

Revised 06/2025

Roche System/Test	LOINC Name	LOINC Number	Component	Property	Time	System	Scale	Method
<b>cobas</b> ® SARS-CoV-2 & Influenza A/B Nucleic acid test for use on the <b>cobas</b> ® <b>liat</b> System	Influenza virus A and B and SARS-CoV-2 (COVID-19) and SARS-related CoV RNA panel - Respiratory specimen by NAA with probe detection	95380-2	Influenza virus A & Influenza virus B & SARS coronavirus 2 & SARS-related coronavirus RNA panel	n/a	Pt	Respiratory	n/a	Probe.amp.tar
<b>cobas</b> SARS-CoV-2 & Influenza A/B Nucleic acid test for use on the <b>cobas</b> <b>liat</b> System	Influenza virus A RNA [Presence] in Respiratory specimen by NAA with probe detection	92142-9	Influenza virus A RNA	PrThr	Pt	Respiratory	Ord	Probe.amp.tar
<b>cobas</b> SARS-CoV-2 & Influenza A/B Nucleic acid test for use on the <b>cobas</b> <b>liat</b> System	Influenza virus B RNA [Presence] in Respiratory specimen by NAA with probe detection	92141-1	Influenza virus B RNA	PrThr	Pt	Respiratory	Ord	Probe.amp.tar
<b>cobas</b> SARS-CoV-2 Nucleic acid test for use on the <b>cobas</b> <b>liat</b> System	SARS-CoV-2 (COVID-19) RNA [Presence] in Respiratory specimen by NAA with probe detection	94500-6	SARS coronavirus 2 RNA	PrThr	Pt	Respiratory	Ord	Probe.amp.tar
<b>cobas</b> C. Difficile Nucleic acid test for use on the <b>cobas</b> <b>liat</b> System	Clostridoides difficile toxin B tcdB gene [Presence] in Stool by NAA with probe detection	74822-8	Clostridoides difficile toxin B tcdB gene	PrThr	Pt	Stool	Ord	Probe.amp.tar
<b>cobas</b> Influenza A/B & RSV Nucleic acid test for use on the <b>cobas</b> <b>liat</b> System	Influenza virus A and B and Respiratory syncytial virus RNA panel - Upper respiratory specimen by NAA with probe detection	85476-0	Influenza virus A & B & Respiratory syncytial virus RNA panel	n/a	Pt	Respiratory. upper	n/a	Probe.amp.tar
<b>cobas</b> Influenza A/B & RSV Nucleic acid test for use on the <b>cobas</b> <b>liat</b> System	Respiratory syncytial virus RNA [Presence] in Respiratory specimen by NAA with probe detection	92131-2	Respiratory syncytial virus RNA	PrThr	Pt	Respiratory	Ord	Probe.amp.tar
<b>cobas</b> Group A Streptococcus Nucleic acid test for use on the <b>cobas</b> <b>liat</b> System	Streptococcus pyogenes DNA [Presence] in Throat by NAA with probe detection	60489-2	Streptococcus pyogenes DNA	PrThr	Pt	Thrt	Ord	Probe.amp.tar
<b>cobas</b> <b>liat</b> SARS-CoV-2, Influenza A/B & RSV nucleic acid test	Influenza virus A and B and SARS-CoV-2 (COVID-19) and Respiratory syncytial virus RNA panel - Respiratory system specimen by NAA with probe detection	95941-1	Influenza virus A & Influenza virus B & SARS coronavirus 2 & Respiratory syncytial virus RNA panel	N/A	Pt	Respiratory	N/A	Probe.amp.tar
<b>cobas</b> <b>liat</b> CT/NG nucleic acid test	Chlamydia trachomatis and Neisseria gonorrhoeae rRNA panel - Urine by NAA with probe detection	70161-5	Chlamydia trachomatis & Neisseria gonorrhoeae rRNA panel	N/A	Pt	Urine	N/A	Probe.amp.tar
<b>cobas</b> <b>liat</b> CT/NG nucleic acid test	Chlamydia trachomatis and Neisseria gonorrhoeae rRNA panel - Specimen by NAA with probe detection	64017-7	Chlamydia trachomatis & Neisseria gonorrhoeae rRNA panel	N/A	Pt	XXX	N/A	Probe.amp.tar
<b>cobas</b> <b>liat</b> CT/NG and/or CT/NG/MG nucleic acid test	Chlamydia trachomatis rRNA [Presence] in Urine by NAA with probe detection	42931-6	Chlamydia trachomatis rRNA	N/A	Pt	Urine	N/A	Probe.amp.tar
<b>cobas</b> <b>liat</b> CT/NG and/or CT/NG/MG nucleic acid test	Chlamydia trachomatis rRNA [Presence] in Specimen by NAA with probe detection	43304-5	Chlamydia trachomatis rRNA	N/A	Pt	XXX	N/A	Probe.amp.tar
<b>cobas</b> <b>liat</b> CT/NG and/or CT/NG/MG nucleic acid test	Neisseria gonorrhoeae rRNA [Presence] in Urine by NAA with probe detection	60256-5	Neisseria gonorrhoeae rRNA	N/A	Pt	Urine	N/A	Probe.amp.tar
<b>cobas</b> <b>liat</b> CT/NG and/or CT/NG/MG nucleic acid test	Neisseria gonorrhoeae rRNA [Presence] in Specimen by NAA with probe detection	43305-2	Neisseria gonorrhoeae rRNA	N/A	Pt	XXX	N/A	Probe.amp.tar
<b>cobas</b> <b>liat</b> CT/NG/MG nucleic acid test	Mycoplasma genitalium rRNA [Presence] in Urine by NAA with probe detection	100703-8	Mycoplasma genitalium rRNA	N/A	Pt	Urine	N/A	Probe.amp.tar
<b>cobas</b> <b>liat</b> CT/NG/MG nucleic acid test	Mycoplasma genitalium rRNA [Presence] in Specimen by NAA with probe detection	100706-1	Mycoplasma genitalium rRNA	N/A	Pt	XXX	N/A	Probe.amp.tar