

cobas[®] Liat[®] Transport Medium Overview



Collection kits direct from Vendor

To use with the **cobas®** Influenza A/B & Influenza A/B & RSV

Supplier	Product name	P/N	Units	Includes
BD	Universal Viral Transport 3 mL Collection Kit with Flexible <u>Minitip</u> Flocked Swab	220531	50	- UTM media vial - Swab
BD	Universal Viral Transport vial, 3mL	220220	50	- UTM/UVT media vial
BD	Flexible <u>Minitip</u> Flocked Swab	220252	100	- Swab
Copan	Copan UTM Tube and Flocked Swab	305C	6 x 50	- UTM media vial - Swab
Copan	Copan UTM – Bulk Tubes (no swabs)	330C	6 x 50	- UTM media vial
Copan	Flexible Nasopharyngeal flocked swab, Nylon tip	503cs01	10 x 50	- Swab

To use with the **cobas®** Strep A

Supplier	Product name	P/N	Units	Includes
BD	BD <u>ESwab</u> Collection Kit Regular Flocked Swab	220245	50	- Liquid <u>Amies</u> vial - Swab
Copan	COPAN <u>ESWAB™</u>	480C	10 x 50	- Liquid <u>Amies</u> vial - Swab



Product Links for Flu | BD

<http://www.bd.com/ds/productCenter/220531.asp>

<http://www.bd.com/ds/productCenter/220220.asp>

<http://www.bd.com/ds/productCenter/220252.asp>

http://www.bd.com/ds/technicalCenter/promotionalIFlyers/ss_1_222804_uvt.pdf



Product Links for Flu | Copan

<http://www.copanusa.com/products/collection-transport/utm-viral-transport/>

UTM™ Tubes in Bulk

Catalog Number: 330C

Tube Size: 16x100mm

Product Description: Conical tube filled with 3mL UTM™ medium

Unit of Measure (UOM): 50 tubes per shelf pack - 6 x 50 tubes per box



Catalog Number: 305C

Tube Size: 16x100mm

Swab Type: Flocked

Product Description: Conical tube filled with 3mL UTM™ medium packaged with one flexible minitip FLOQSwab™, sterile

Unit of Measure (UOM): 50 kits per shelf pack - 6 x 50 kits per box



http://www.copanusa.com/files/1814/2352/3834/Copan_FLOQSwab_broch-web.pdf

Flexible Mini Tip



Product Links for Strep | BD and Copan

<http://www.copanusa.com/products/col-lection-transport/eswab/>



<http://www.bd.com/ds/productCenter/20245.asp>

