


DISCOVERY Detection Reagents for IHC

DISCOVERY detection reagents deliver the flexibility and quality that your research demands. Our menu includes biotin-free detection such as DISCOVERY HQ HRP and the DISCOVERY NP AP detection, featuring our unique hapten-linked multimer technology to maximize your signal to noise ratio, and a variety of distinct chromogens for multiplexing. We build the tools you need to piece together answers to the toughest questions in research.




HRP multimer detection systems			
Product name	Catalog number	Ordering code	Number of tests
DISCOVERY HQ HRP hapten-linked multimer detection			
Multimer only			
■ DISCOVERY Anti-HQ HRP	760-4820	07017936001	125 tests
■ DISCOVERY Anti-Mouse HQ	760-4814	07017782001	125 tests
■ DISCOVERY Anti-Rabbit HQ	760-4815	07017812001	125 tests
OmniMap HRP multimer detection			
Multimer only			
■ DISCOVERY OmniMap anti-Gt HRP (RUO)	760-4647	06607233001	50 tests
■ DISCOVERY OmniMap anti-Ms HRP (RUO)	760-4310	05269652001	125 tests
■ DISCOVERY OmniMap anti-Rb HRP (RUO)	760-4311	05269679001	125 tests
■ DISCOVERY OmniMap anti-Rt HRP (RUO)	760-4457	05891892001	125 tests
Chromogen and multimer			
■ ■ OmniMap DAB anti-Ms detection kit	760-150	05266556001	125 tests
■ ■ OmniMap DAB anti-Rb detection kit	760-149	05266548001	125 tests
■ ■ OmniMap DAB anti-Gt (RUO)	760-182	06525725001	125 tests
■ ■ OmniMap DAB anti-Rt (RUO)	760-158	06525709001	125 tests
UltraMap HRP multimer detection			
Multimer only			
■ DISCOVERY UltraMap anti-Gt HRP (RUO)	760-4648	06607241001	50 tests
■ DISCOVERY UltraMap anti-Ms HRP (RUO)	760-4313	05269695001	125 tests
■ DISCOVERY UltraMap anti-Rb HRP (RUO)	760-4315	05269717001	125 tests
■ DISCOVERY UltraMap anti-Rt HRP (RUO)	760-4456	05891884001	125 tests
Chromogen and multimer			
■ ■ UltraMap DAB anti-Ms detection kit	760-152	05266572001	125 tests
■ ■ UltraMap DAB anti-Rb detection kit	760-151	05266564001	125 tests
■ ■ UltraMap DAB anti-Rt (RUO)	760-157	06525717001	125 tests
■ ■ UltraMap DAB anti-Gt (RUO)	760-183	06525733001	
HRP-driven chromogens			
■ DISCOVERY ChromoMap DAB kit (RUO)	760-159	05266645001	125 tests
■ DISCOVERY Purple kit (RUO) *	760-229	07053983001	50 tests
■ DISCOVERY Silver kit (RUO) *	760-227	07053649001	50 tests
DISCOVERY Inhibitor *	760-4840	07017944001	250 tests
■ DISCOVERY Teal HRP Kit*	760-247	08254338001	50 tests

* DISCOVERY Inhibitor is required when using DISCOVERY Purple and DISCOVERY Silver chromogens




Any of these biotin-free HRP multimer detection kits can be used with any of these chromogens

Amplification			
Product name	Catalog number	Ordering code	Number of tests
DISCOVERY Amp HQ kit	760-052	06472320001	50 tests
DISCOVERY Amplification anti-HQ HRP multimer	760-4602	06442544001	50 tests
DISCOVERY Amplification anti-HQ AP multimer	760-4521	06472311001	50 tests



Amplification in any color can be done starting with an HRP multimer detection kit

AP multimer detection systems			
Product name	Catalog number	Ordering code	Number of tests
DISCOVERY NP AP hapten-linked multimer detection			
DISCOVERY Anti-Mouse NP	760-4816	07425309001	125 tests
DISCOVERY Anti-Rabbit NP	760-4817	7425317001	125 tests
DISCOVERY Anti-NP AP	760-4827	7425325001	125 tests
UltraMap AP multimer detection			
Multimer only			
DISCOVERY UltraMap anti-Ms Alk Phos	760-4312	5269687001	125 tests
DISCOVERY UltraMap anti-Rb Alk Phos	760-4314	5269709001	125 tests
Chromogen and Multimer			
UltraMap Red anti-Ms detection kit	760-154	5266602001	125 tests
UltraMap Red anti-Rb detection kit	760-153	5266599001	125 tests
UltraMap Blue anti-Ms detection kit	760-156	5266637001	125 tests
UltraMap Blue anti-Rb detection kit	760-155	5266629001	125 tests
AP-driven chromogens			
DISCOVERY Red Kit (RUO)	760-228	07425333001	125 tests
DISCOVERY ChromoMap Blue kit (RUO)	760-161	5266661001	125 tests
DISCOVERY Yellow kit (RUO)	760-239	07698445001	125 tests



Any of these biotin-free AP multimer detection kits can be used with any of these chromogens

Biotin-based detection systems			
Product name	Catalog number	Ordering code	Number of tests
DISCOVERY DAB Map detection kit (RUO)	760-124	05266360001	250 tests
DISCOVERY RedMap kit (RUO)	760-123	05266351001	250 tests
DISCOVERY BlueMap kit (RUO)	760-120	05266327001	250 tests

DISCOVERY products are for research use only. Not for use in diagnostic procedures.

VENTANA and DISCOVERY are trademarks of Roche.
All other product names and trademarks are the property of their respective owners.

© 2018 Roche MC-US-00677-0718

Roche Diagnostics Corporation
9115 Hague Road
Indianapolis, Indiana 46256

usdiagnostics.roche.com