









DISCOVERY Detection Reagents for ISH

DISCOVERY Detection reagents deliver the flexibility and quality that your research demands. Our menu includes our biotin-free ISH specific DISCOVERY anti-DIG and anti-DNP Multimers designed for full automation. These multimers may be used with any of our chromogenic or fluorescent detection kits to generate the signal required for your assay. We build the tools you need to piece together answers to the toughest questions in research.



HRP multimer detection for DIG and DNP labeled probes

Product name	Catalog number	Ordering code	Number of tests
HRP Multimer detection			
 DISCOVERY anti-DIG HRP Multimer (RUO)	760-4822	07256299001	50 tests
 DISCOVERY anti-DNP HRP Multimer (RUO)	760-4821	07339712001	50 tests
HRP driven chromogens			
 DISCOVERY ChromoMap DAB Kit (RUO)	760-159	05266645001	125 tests
 DISCOVERY Purple Kit (RUO)*	760-229	07053983001	50 tests
 DISCOVERY Silver Kit (RUO)*	760-227	07053649001	50 tests
DISCOVERY Inhibitor*	760-4840	07017944001	250 tests
HRP driven fluorophores			
 DISCOVERY Rhodamine Kit (RUO)*	760-233	07259883001	50 tests
 DISCOVERY FITC Kit (RUO)*	760-232	07259212001	50 tests
DISCOVERY Inhibitor*	760-4840	07017944001	250 tests



Any of these Biotin-free HRP multimer detection kits for ISH can be used with any of these chromogens or fluorescence detection kits


*DISCOVERY Inhibitor is required when using DISCOVERY Purple, DISCOVERY Silver, DISCOVERY FITC, or DISCOVERY Rhodamine Kits

AP multimer detection for DIG and DNP labeled probes


Product name	Catalog number	Ordering code	Number of tests
DISCOVERY anti-DIG AP Multimer (RUO)	760-4825	07256302001	50 tests
DISCOVERY anti-DNP AP Multimer (RUO)	760-4824	07339739001	50 tests

AP driven Chromogens

DISCOVERY ChromoMap Red Kit (RUO)	760-160	05266653001	125 tests
DISCOVERY ChromoMap Blue Kit (RUO)	760-161	05266661001	125 tests



Any of these Biotin-free AP multimer detection kits can be used with any of these chromogens



ROS 5' (cent) DIG Probe (RUO) with DISCOVERY Rhodamine Kit and ROS 3' (telo) DNP Probe (RUO) with DISCOVERY FITC Kit in tonsil

Alternatively a secondary linking antibody may be used along with any of our Biotin-free anti species multimer detection kits.

HRP multimer detection systems			
Product name	Catalog number	Ordering code	Number of tests
DISCOVERY HQ HRP multimer detection			
Multimer only			
DISCOVERY Anti-HQ HRP	760-4820	07017936001	125 tests
DISCOVERY Anti-Mouse HQ	760-4814	07017782001	125 tests
DISCOVERY Anti-Rabbit HQ	760-4815	07017812001	125 tests

UltraMap HRP multimer detection			
Multimer only			
DISCOVERY UltraMap anti-Gt HRP (RUO)	760-4648	06607241001	50 tests
DISCOVERY UltraMap anti-Ms HRP (RUO)	760-4313	05269695001	125 tests
DISCOVERY UltraMap anti-Rb HRP (RUO)	760-4315	05269717001	125 tests
DISCOVERY UltraMap anti-Rt HRP (RUO)	760-4456	05891884001	125 tests

OmniMap HRP multimer detection			
Multimer only			
DISCOVERY OmniMap anti-Gt HRP (RUO)	760-4647	06607233001	50 tests
DISCOVERY OmniMap anti-Ms HRP (RUO)	760-4310	05269652001	125 tests
DISCOVERY OmniMap anti-Rb HRP (RUO)	760-4311	05269679001	125 tests
DISCOVERY OmniMap anti-Rt HRP (RUO)	760-4457	05891892001	125 tests

Secondary linking antibody			
Rabbit Anti-DNP Antibody	780-4335	05273447001	100 tests

AP multimer Detection systems			
Product name	Catalog number	Ordering code	Number of tests
UltraMap AP multimer detection			
Multimer only			
DISCOVERY UltraMap anti-Ms Alk Phos	760-4312	05269687001	125 tests
DISCOVERY UltraMap anti-Rb Alk Phos	760-4314	05269709001	125 tests

DISCOVERY products are for research use only. Not for use in diagnostic procedures.

© 2018 Roche MC-US-01567-0918

Roche Diagnostics Corporation
9115 Hague Road
Indianapolis, Indiana 46256

usdiagnostics.roche.com

VENTANA and DISCOVERY are trademarks of Roche.
All other product names and trademarks are the property of their respective owners.