



DISCOVERY detection reagents for IHC



Multimer detection HRP/AP kits

| MULTIMER No need for biotin blocking | | | |
|--|--|-----------------|-----------|
| OmniMAP | | UltraMAP | |
| HRP | | HRP | AP |
| Ms, Rb, Rt, Gt | | Ms, Rb, Rt, Gt | Ms, Rb |




AP
Chromogens:
 Red, Yellow, Blue
No fluorophores

HRP
Chromogens:
 Purple, Silver, Teal,
 DAB, Green
Fluorophores: DCC, Cy5,
 FAM, Red610, Rhodamine,
 FITC, Rhodamine 6G

| Multimer detection HRP/AP kits | | | | |
|---------------------------------------|---------|----------|-------------|-----------|
| Multimer HRP | | | | |
| OmniMAP HRP | Anti-Ms | 760-4310 | 05269652001 | 125 tests |
| | Anti-Rb | 760-4311 | 05269679001 | 125 tests |
| | Anti-Gt | 760-4647 | 06607233001 | 50 tests |
| | Anti-Rt | 760-4457 | 05891892001 | 125 tests |
| UltraMAP HRP | Anti-Ms | 760-4313 | 05269695001 | 125 tests |
| | Anti-Rb | 760-4315 | 05269717001 | 125 tests |
| | Anti-Gt | 760-4648 | 06607241001 | 50 tests |
| | Anti-Rt | 760-4456 | 05891884001 | 125 tests |
| Multimer AP | | | | |
| UltraMAP AP | Anti-Ms | 760-4312 | 05269687001 | 125 tests |
| | Anti-Rb | 760-4314 | 05269709001 | 125 tests |

Hapten detection kits

| HAPTEN | |
|---|---------------|
| Increased sensitivity, decreased background | |
| NP | HQ |
| AP Ms, Rb | HRP Ms, Rb |




| Hapten detection kits | | | | |
|-----------------------|-----------------------|----------|-------------|-----------|
| HQ Detection | DISCOVERY Anti-Ms HQ | 760-4814 | 07017782001 | 125 tests |
| | DISCOVERY Anti-Rb HQ | 760-4815 | 07017812001 | 125 tests |
| | DISCOVERY Anti-HQ HRP | 760-4820 | 07017936001 | 125 tests |
| NP Detection | DISCOVERY Anti-Ms NP | 760-4816 | 07425309001 | 125 tests |
| | DISCOVERY Anti-Rb NP | 760-4817 | 07425317001 | 125 tests |
| | DISCOVERY Anti-NP AP | 760-4827 | 07425325001 | 125 tests |

HQ/BF amplification kits

DISCOVERY HQ/BF amplification

Increase signal intensity and reduce background



DISCOVERY amplification requires a HRP multimer to drive the tyramide reaction

HQ/BF amplification kits

| | | | | |
|-------------------------|---------------------------|----------|-------------|----------|
| HQ Amplification | DISCOVERY HQ Amp kit | 760-052 | 06472320001 | 50 tests |
| | DISCOVERY Amp Anti-HQ AP | 760-4521 | 06472311001 | 50 tests |
| | DISCOVERY Amp Anti-HQ HRP | 760-4602 | 06442544001 | 50 tests |
| BF Amplification | DISCOVERY BF Amp kit | 760-226 | 07529422001 | 50 tests |
| | DISCOVERY Amp Anti-BF HRP | 760-4825 | 07529449001 | 50 tests |

Chromogen and Fluorophore kits

Chromogen kits

| | | | |
|--------------------------------|---------|-------------|-----------|
| DISCOVERY Red kit | 760-228 | 07425333001 | 50 tests |
| DISCOVERY Yellow kit | 760-239 | 07698445001 | 50 tests |
| ChromoMAP Blue kit | 760-161 | 05266661001 | 50 tests |
| ChromoMAP Red kit | 760-160 | 05266653001 | 50 tests |
| ChromoMAP DAB kit | 760-159 | 05266645001 | 125 tests |
| DISCOVERY Purple kit | 760-229 | 07053983001 | 50 tests |
| DISCOVERY Silver kit | 760-227 | 07053649001 | 50 tests |
| DISCOVERY Teal HRP kit | 760-247 | 08254338001 | 50 tests |
| DISCOVERY Green HRP kit | 760-271 | 08478295001 | 50 tests |

Fluorophore kits

| | | | |
|-----------------------------------|---------|-------------|----------|
| DISCOVERY DCC kit | 760-240 | 07988192001 | 50 tests |
| DISCOVERY FAM kit | 760-243 | 07988150001 | 50 tests |
| DISCOVERY Red610 kit | 760-245 | 07988176001 | 50 tests |
| DISCOVERY Cy5 kit | 760-238 | 07551215001 | 50 tests |
| DISCOVERY Rhodamine 6G kit | 760-244 | 07988168001 | 50 tests |
| DISCOVERY Rhodamine kit | 760-233 | 07259883001 | 50 tests |
| DISCOVERY FITC kit | 760-232 | 07259212001 | 50 tests |

DISCOVERY products are for research use only. Not for use in diagnostic procedures.

VENTANA and DISCOVERY are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

© 2020 Roche. MC-US-06771-0520

Roche Diagnostics Corporation
9115 Hague Road
Indianapolis, Indiana 46256

diagnostics.roche.com